





Instructions for use Models S 280, S 380



# Help

### **Hyperlinks**

Hyperlinks are cross-references in the text. They are marked with the symbol. Entries in the table of contents also have hyperlink functionality.

- Clicking takes you to the relevant page.
- Clicking in the Acrobat navigation bar takes you back to the start page.

#### Full-text search for a term

Clicking pens the "Search" window.

You can search for any term in this Acrobat document or optionally search all of the Acrobat files stored on your computer or those that are accessible.

#### Other interactive functions

The Adobe Reader and all Adobe Acrobat program versions offer the experienced user a whole raft of highly useful functions for interactive working with documents. For further information see the corresponding user manuals or online help.

Display table of contents

Jump to Help page (this page)

Contents 2 of 88 Help

# Version, Imprint, Type label

Instructions for use G77 Version 1.0 (06.07.2009) for Völker care bed models S 280, S 380 built after April 2005

© by Völker AG 2009

Völker AG Wullener Feld 79 58454 Witten

#### **GERMANY**

Tel.: +49 2302 96096-0 Fax: +49 2302 96096-16 E-mail: info@voelker.de

All rights reserved. Reproduction - even in part - is not permitted


in part - is not permitted.

We reserve the right to make changes to reflect technical advances.

The contents of this document are subject to change without prior notification.

Customers are advised to contact the responsible area sales manager before placing an order.

### Type label



The type label is located on the inside of the head panel.

Raise the rear section to read the type label.



### Table of contents

Activation 27 Help 🕮 2 Cleaning and disinfection 49 Version, Imprint, Type label 2 3 General operating instructions 28 Cleaning 50 Table of contents 4 Preparation 29 Wipe and spray disinfection \$\omega\$ 50 Foreword 6 Electrical activation A 30 Maintenance 4 52 Using the battery pack (option) 4 31 Notes P 7 Taking out of service 4 32 Staff training 4 53 General information A 8 Functional check A 33 Safety notes 4 53 Designated purpose 49 Maintenance schedule 4 54 General regulations and Operation 34 user training / instruction Key lock 235 Technical check \$\Pi\$ 55 **10** Central castor adjustment 4 36 Visual inspection 4 56 Side rails 2 37 General safety notes 4 11 Functional check of the side rails General safety notes 2 37 **56** Functional description 15 Model S 280 🕮 38 Functional check of the brakes 4 57 Overview 16 Model S 380 🕮 39 Functional check of the motors Model S 280 <sup>□</sup> 16 Back section 41 **57** Model S 380 P 17 Mechanical rapid lowering of the back Mains connection cable 4 57 Hand control with hook 118 section / CPR function (option) 42 Cabling 57 Upper leg section 43 Housing 4 57 Transverse hand control (option for S 380) 🕮 19 Lower leg section 44 Mechanical check 4 57 Nurse keypad with complete lock Lying surface height 45 Changing the battery 4 57 (option) Decided 20 Anti-Trendelenburg and Trendelenburg Measurement as per VDE 0751-1 Trapeze bar and accessory holders, positioning 46 **57** accessory rail 22 Comfort seating position 47 Checking the trapeze bar grab Versions and options 23 Bed extension (option) 48 handle III 58 Accessories 25

Contents 4 of 88 Help

Other accessories 58

Troubleshooting 59

Table of faults 60
Service points 62

Appendix 63
Symbols used 64
Technical data 65
Classification 67

Service life / disposal 🕮 70

Manufacturer's declarations, forms, electronic instructions for use □ 71

Type label 48

Contents 5 of 88 Help

### **Foreword**

We are delighted that you have chosen Völker care beds.

We are also grateful for the trust you have placed in our company and our products.

This step is undoubtedly the result of extensive considerations and examinations of the requirements that you want to make on new care beds based on your previous experiences.

You clearly had good reasons to choose Völker care beds.

We promise you that Völker care beds will not disappoint you.

Not for nothing do Völker care beds enjoy a worldwide reputation as particularly innovative medical aids. This refers not only to the design principle, which was completely redeveloped by Völker. It also refers to the multitude of product advantages which were continually tested for their practicability in practice and improved. These will support the comfort of the occupant/patient and also help to reduce the load of daily care work.

Now, every care bed has product features that are of practical benefit to their users. However, as far as we know, none of them offers the range of advantages that a Völker care bed does.

Völker care beds not only look great, but they also offer functions that can be controlled or adjusted mechanically, although most of them are controlled or adjusted using electric motors or electronic components.

If you purchase these beds, the responsibility for their correct and appropriate use is transferred to you. Consequently, we strongly advise that you consult the enclosed instructions for use to learn about the technical features, handling and use of all the functions.

We wish you every success with Völker care beds

Heinrich Völker Chairman, Völker Aktiengesellschaft

### **Notes**



# **Notes | General information**

#### **General notes**

You have purchased a care bed from Völker AG. This care bed has been built in accordance with the applicable national and international standards and the regulations reflecting the current state of technology.

Völker care beds satisfy all the requirements in terms of safety and functionality. They are tested in accordance with international standards and bear the CE mark, which documents the beds' compliance with essential requirements for medical products.

Please read the general safety notes

11. Please also note (with particular attention to any warranty claims) the further notes on the following pages.

### Standard design

The standard design of the bed can be supplied with various options. A description of these options can be found in the section entitled Versions and options 23.

### **Copyright protection**

These instructions for use may only be transferred to third parties with the written consent of Völker AG. All documents are protected under copyright.

### Warranty and liability

Völker AG is liable for any faults or failures, not including further claims arising in the context of the warranty obligations detailed in the main agreement. Claims for compensation, for whatever legal reason such claims may be raised, will not be entertained.

We reserve the right to make technical modifications as part of the further development of the care beds which form the subject of these instructions for use.

We accept no liability for damage and operational faults caused as a result of misuse and/or non-observance of these instructions for use.

The portrayal of the accessories does not necessarily match the actual product.

# **Notes | Designated purpose**

#### Correct use

Völker care bed models S 280 and S 380 are designed for the laying down and care of occupants/patients in care institutions, hospitals and in suitable rooms in residential buildings.

The bed is intended for use by people over the age of 12 or who are taller than 146 cm.

The safe working load of the bed is 230 kg (210 kg for beds built before 10/2008. The maximum permissible weight of the individual is derived by deducting the weight of the mattress, the trapeze handle and other accessories from this.

Any use of the Völker care bed other than for the purpose intended excludes the company from any possible liability.

### Inappropriate use

Inappropriate use can be dangerous. This includes, but is not limited to:

- Incorrect actuation of electrical functions and uncontrolled positioning;
- Operation of the care bed by the occupant/patient without having received prior instruction in how to do so;
- Use of electrical equipment on the bed that is not intended for such use (subject to the operator's obligation to exercise care);
- · Pulling on cables to move the bed;
- Loosening electrical plug connections by pulling on the cable;
- Use of the bed on a slope of more than ten degrees of inclination (the bed's brakes are designed for an angle of inclination of no more than ten degrees);
- Any attempt to move the bed while it is in braked position;

- Use of the bed for transport with a vehicle;
- Overloading of the bed beyond the specified safe working load.



**CAUTION** If, in an emergency situation, it is impossible to avoid putting children under the age of 12 or people who are less than 146 cm tall in the bed, protective covers must be placed on the side rails. This also applies to the use of the bed by weak or confused patients. The use of the bed for children under the age of 8 is absolutely forbidden.

# Notes | General regulations and user training / instruction

### **General regulations**

The care bed must only be operated and used in accordance with its designated purpose, in line with the conditions of the Medical Products Directive (MPG) and approved legislation pursuant to this, the generally-acknowledged rules of technology and the stipulations of occupational safety and accident protection guidelines. The care bed must not be operated in a faulty state that could endanger its occupant/patient, care personnel or third parties.

#### **User training**

The care bed may only be operated by individuals whose training or understanding and experience offer surety for correct handling (medical products directive).

#### **User instruction**

The thorough introduction of care personnel in the operation of the bed can be provided by Völker or its representative at the customer's request.

Attendance of such training can be certified and confirmed by Völker using the form provided for this purpose, specifying the name, date and signature.

Occupants/patients must be instructed in the use of the bed before care personnel hand over the hand control to them.

#### Other requirements

Whoever is in charge of the activation, operation or preparation of the bed must have been given a copy of these instructions for use (in printed or electronic form) and have read them.

To avoid operating errors and to safeguard the smooth operation of the bed, care personnel must always have access to the safety notes below.

# Notes | General safety notes 1/4



Warning symbols Information marked with this symbol must be read and its content strictly observed.



**DANGER** represents an immediate threat of danger that can cause serious physical injury or death.



WARNING represents potentially dangerous situations that can lead to serious physical injury or death.



**CAUTION** represents potentially dangerous situations that can cause slight physical injuries.

**NOTE** warns of potential damage to objects or property.

#### Before first activation

Before the care bed is put into action for the first time, care personnel must read these instructions for use in full and with care.

Before the bed is activated for the first time, care personnel must be instructed in the handling of the bed using the instructions for use. The potential dangers that can arise despite correct operation of the bed must also be pointed out in full

### Before and during use

Before each use of the bed, the user must be sure that the bed is in a good, safe condition and that safe use is ensured (Functional check (4) 33).

#### Position of the care bed



**CAUTION** To avoid injuries caused by falling, we recommend (except while care is being given) that the bed is generally set at its lowest position with the brake activated.

### Transporting the bed



**CAUTION** When transporting the bed, it must always be ensured that the mains connection cable is not touching the floor. The bed should only be moved on a solid floor. Never attempt to push it over obstacles of over 2 cm in height. The maximum angle of inclination of the floor must not exceed 10°.

# Notes | General safety notes 2/4

#### Four castor central braking



#### **CAUTION "Risk of accident"**

If the bed is not being transported, the castors must always be in the braked position, since the bed may be required as a support for when the occupant/patient stands up or lies down. If the bed rolls away while unbraked, this can lead to a serious fall. After activating the central holding brake, as 36 it must be checked that the bed is actually fixed, i.e. the castors are adequately braked.

The bed can be in a non-fullybraked status even after first activation or reactivation, and consequently it must be checked that the castors are correctly braked.

#### One-sided load on the bed



**NOTE** In order to prevent onesided loads on the bed, it must not be used as a seat for persons other than the occupant/patient (i.e. visitors must not sit on the edge of the bed).

#### Side rails



### WARNING "Risk of entrapment" In the case of occupants whose physical or mental condition makes necessary to use side rails

se physical or mental condition makes necessary to use side rails to protect them from falling out of bed, the following safety measures must be observed:

- The legal permissibility of using side rails must be ascertained.
- The side rails may only be operated by trained care personnel.

- Ensure that the side rails (or parts thereof) are either fully raised and locked in position or completely lowered.
- It must be ensured that the occupant does not come into contact with the side rail elements when the electrical lying surface adjustment mechanism is actuated. It is also important to ensure that no part of the body is sticking out through the side rails.
- If the side rails are used with a person whose psychological condition makes their use necessary, then it must be ensured that the hand control is kept out of their reach or its functions are locked. It is strongly recommended that side rail covers are used.

# Notes | General safety notes 3/4

 If these safety measures are not observed by care personnel, injuries can be caused to hands, knees, fingers, feet, legs and hips, along with haematomas and other injuries as a result of entrapment. In case of children or people who are less than 146 cm tall, nonobservance of these guidelines can lead to death!



WARNING "Risk of entrapment" On model S 280, the side rails must either be fully raised and securely locked in position, or be completely lowered.



WARNING "Danger of injury"
If the side rails are damaged, the bed must not be used and must be repaired.

### Height adjustment



DANGER "Risk of entrapment between the lower frame and/or floor and the bed frame when the bed is lowered"

It must be ensured that no people, limbs, pets, bed linen or other objects are caught between the bed frame and the lower frame and/or floor.



### DANGER "Danger of movement"

If any movement of the bed could represent a danger to the occupant/patient, all functions must be locked.

#### **Accessories**



### WARNING "Risk of injury"

Only original Völker accessories should be used! Third-party accessories must be subjected to testing before use.

### **Using lifting devices**



### WARNING "Risk of injury"

No lifting device must be fastened directly to the bed (patient transport, repair).

The lifting devices specified are appliances that can be attached to the bed for transport purposes. Patient lifters can be used.

# Notes | General safety notes 4/4

#### Using oxygen equipment



#### **DANGER "Risk of fire"**

Do not use any oxygen equipment other than that which is administered via nasal prongs or masks. Do not use the bed in a room with a potential explosion risk. (Provided it has been excluded (e.g. based on information in the instructions for use of the equipment being used) that the use of the equipment may cause the O<sub>2</sub> concentration to rise to such a degree that there is an explosion risk (even when there is no fault), then the equipment can be used).

#### Rail spacers

When using the rail spacer, please read the separate instructions for use for this accessory. During technical checks, the rail spacers should also be checked to ensure they are suitable for the size of side rail used.

#### Cleaning and disinfection

In order to maintain consistent functioning, the care bed should be cleaned, disinfected and tested as soon as possible following each use, so that it can be reused immediately without risk.

Dangers can arise from the incorrect cleaning / disinfection  $\square$  48 of the bed.

### Maintenance and repair

Anyone responsible for carrying out maintenance and repair work must at least have read the safety notes and the service manual and be qualified in accordance with MPBetreibV §§ 4 and 6.

After maintenance work or repairs have been carried out, a technical check 55 must be carried out on the affected parts and/or functions. During this check, it must be determined that the bed can be used in accordance with the specifications without risk to the occupant/patient, user or a third party.

The technical check must be carried

out at least once a year and after every lengthy period of non-use.

Any discernible damage, such as signs of wear and tear, loose screws or breaks/fractures must be eliminated immediately.

### Electromagnetic / static interference

The care beds in model series S 280 and S 380 satisfy the EMC requirements in accordance with the law on medical products (MPG). The basis for testing is standard EN 60601-1-2.

# **Functional description**

The features of the Völker care bed and its function are set out in the section **Functional description**.



### **CONTENTS**

Overview 16

Model S 280 16

Model S 380 17

Hand control with hook 4 18

Transverse hand control (option for S 380) 4 19

Nurse keypad with complete block (option) 20

Trapeze bar and accessory holders, accessory rail 22

Versions and options 23


Accessories 25

# Functional description | Overview | Model S 280



Contents 16 of 88 Help

# Functional description | Overview | Model S 380



Contents 17 of 88 Help

# Functional description | Hand control with hook

Back section up

Upper leg section up

Anti-Trendelenburg positioning <sup>1</sup>

Lying surface up



Back section down

Upper leg section down

Auto-Contour <sup>2</sup> (optional)

Lying surface down





Hand control unlocked



Hand control locked

**WARNING** When actuating motorised adjustments with the side rails raised, it must be ensured that the occupant/patient does not have any contact with the side rails, and that no parts of his or her body or that of other people are sticking through the side rails or are trapped between the lying surface and the lower frame and/or floor!

Contents 18 of 88 Help

 $<sup>\</sup>triangle$ 

<sup>&</sup>lt;sup>1</sup>Head end raised

<sup>&</sup>lt;sup>2</sup>Back and upper leg section raised

# Functional description | Transverse hand control (option for S 380)

Lying surface Anti-Trendelen- Upper leg sec- Back section up burg positioning <sup>1</sup> tion up up



Lying surface down

Auto-Contour <sup>2</sup> (optional)

Upper leg sec-

Back section down

#### Reverse:



Hand control unlocked



Hand control locked



**WARNING** When actuating motorised adjustments with the side rails raised, it must be ensured that the occupant/patient does not have any contact with the side rails, and that no parts of his or her body or that of other people are sticking through the side rails or are trapped between the lying surface and the lower frame and/or floor!

Contents 19 of 88 Help

<sup>&</sup>lt;sup>1</sup> Head end raised

<sup>&</sup>lt;sup>2</sup> Back and upper leg section raised


# Functional description | Nurse keypad with complete lock 1/2 (option)

Anti-Trendelenburg positioning <sup>1</sup>

Lying surface up <sup>3</sup>

Back section up

Upper leg section up



Trendelenburg positioning <sup>2</sup>

Lying surface down <sup>3</sup>

Back section down

Upper leg section down

- <sup>1</sup>Head end raised
- <sup>2</sup> Head end lowered
- <sup>3</sup> The nurse keypad can be equipped with an automatic function (double-click) on a country-by-country basis (see next page).



WARNING When actuating motorised adjustments with the side rails raised, it must be ensured that the occupant/patient does not have any contact with the side rails, and that no parts of his or her body or that of other people are sticking through the side rails or are trapped between the lying surface and the lower frame and/or floor!

#### Reverse:



Keypad unlocked



In this position, the keypad and hand control are completely locked.

Contents 20 of 88 Help

# Functional description | Nurse keypad with complete block 2/2 (option)

# Double-click function "Lying surface up/down"

The nurse keypad can be equipped with an automatic function (double-click) for lying surface height adjustment on a country-by-country basis.

With a double-click on the "Lying surface up" or "Lying surface down" button, the lying surface moves to the highest or lowest position, respectively. This function can be stopped at any time by pressing any button.



WARNING "Risk of entrapment" If the double-click function is being used, the care giver must supervise the bed occupant/patient until the adjustment procedure has completed.



DANGER "Risk of entrapment between the lower frame and/or floor and the bed frame when the bed is being lowered" It must be ensured that:

No people, limbs, pets, bed linen or other objects are trapped between the bed frame and lower frame and/or floor during adjustment manoeuvres.

# Functional description | Trapeze bar and accessory holders, accessory rail

On the inside of the head panel are holders for the trapeze bar and accessories.

The trapeze bar and other accessories must be slotted into the holders until they audibly engage.



### WARNING "Risk of injury"

Ensure that the trapeze bar is completely slotted into the holder and securely seated.


Note: the safe working load of the trapeze bar is max. 75 kg.



### WARNING "Risk of falling"

The trapeze bar must not be used by the occupant as an aid for climbing into the bed.

The trapeze bar must never jut out beyond the outer edge of the bed and then be used as a pulling-up aid (e.g. when getting out of a wheelchair).



Accessories holder (e.g. for drip stands)

There is an accessory rail on either side of the bed to accommodate accessories.



Accessory rail (e.g. for urine bottle basket, universal hook, etc.)

# Functional description | Versions and Options 1/2

The standard designs of the care beds can be supplied with various versions and options:

### Version/option

Hand control (versions)

### Description

1. With hook (standard) \*:



2. Transverse on the side rail with clip (option for S 380)\*:



### Version/option

Nurse keypad (options)

#### Description

1. In drawer:



2. In bed linen storage area:



3. Attached to the accessory rail:



<sup>\*</sup>A hand control with Auto-Contour button is optionally available

# Functional description | Versions and Options 2/2

Version/option	Description
Castors	The standard design has 150 mm standard castors. Different types of castors are available as an option. The design and the diameter are variable. This can lead to the variation in the lifting adjustment area of around 20 mm.
Bed extension (option)	The S 380 model can be extended by 20 cm with a telescopic bed extension.
Side rails (versions)	Models S 280 and S 380 K (except design version LP) can be equipped with various side rail versions:  Back/lower leg section:
	Can be pulled out up to 34* cm (standard)
	<ol> <li>Can be pulled out up to 40** cm (version S 380)</li> <li>Can be pulled out up to 38** cm (version S 280)</li> </ol>
	<ul> <li>continuously for model S 280.</li> <li>Measured from the top edge of the side rail to the lying surface (without the mattress).</li> </ul>

These instructions for use cover all of the versions and options listed.

Precise details of the supplied bed designs can be found in the order specifications for your beds. If the original bed specification should no longer be available, please contact Völker Customer Services. Make a note of the Völker serial number (ID No.) on the type label 68 before you call.

# Functional description | Accessories 1/2

In order to obtain the greatest possible degree of flexibility, Völker offers a wide range of easy-to-attach accessories. The care beds are equipped as standard with holder devices for accessories, such as drip stands and trapeze bars. Urine bottle baskets, universal hooks, standard bars, etc., can be mounted on the accessory rails provided on both sides of the bed.

Further information about accessories can be found in our current information brochure or on the Internet at www.voelker.de . Our staff will gladly provide you with more details on the accessories that are able for your bed model.

# $\wedge$

**WARNING** Only original Völker accessories should be used! Third-party accessories must be subjected to testing before use.

#### **Mattresses**

Mattress size	Mattress frame size	Density
88 x 200/210/220 x 12 cm	90 x 200/210/220 cm	40 - 50 kg/m <sup>3</sup>
98 x 200/210/220 x 12 cm	100 x 200/210/220 cm	40 - 50 kg/m <sup>3</sup>
108 x 200/210/220 x 12 cm *	110 x 200/210/220 cm *	40 - 50 kg/m <sup>3</sup> *

\* only possible for S 380 model




**DANGER** The height of the raised side rail above the mattress must always be greater than or equal to 220 mm; otherwise the occupant/patient may accidentally fall out of bed. Please note that the height of the mattress has a direct influence on this

To minimise the risk of injury, only use mattresses with the dimensions and properties detailed above. If you do not use Völker mattresses, please contact a dealer in whom you have confidence.



**DANGER** If mattresses are used that do not match these specifications, there is a risk of suffocation!



# Functional description | Accessories 2/2

### Use of securing systems

Securing systems such as belts or straps should only be used exactly as specified by the manufacturer.

If securing systems in the form of abdominal belts are used, then it must be ensured that the side rails are completely raised. The gap in the middle of the rails must always be closed using a rail spacer.



**DANGER** When using securing systems and rail spacers, please note the separate instructions for use pertaining to these accessories.

The lying surfaces must **never** be adjusted while the patient is secured **and** must always be in the lowest position!

The lying surface adjustment functions must be locked when a patient is secured, and the hand control must be kept out of the occupant's/patient's reach!

### **Activation**



Contents 27 of 88 Help

# **Activation | General operating instructions**

#### On-time

The maximum on-time for the electromotive bed functions is specified on the bed (type label 

68) or in the technical data sheet.

1 min/10 min. means that each electromotive adjustment may be operated for a maximum of 1 minute in 10 minutes (protection against overheating).

NOTE Should the maximum on-time of 1 minute be exceeded repeatedly or for longer, safety cutout devices on the bed may cause the electromechanical motor system to fail. The bed must not be manoeuvred using the motors until it has cooled down sufficiently! In the event of severe overloading, damage can occur to the motor.

#### **Battery pack (option)**

The battery pack in the bed has a charge capacity that is equivalent in theory to a constant operation of 2 lifting and lying surface adjustments with a working load of 230 kg.

**NOTE** If the bed is parked at its location and the mains plug is not connected, this will cause self-discharging of the battery pack!

Deep-discharged battery packs can be damaged to the extent that premature replacement is needed!

Appropriate and correct use of the battery is essential for achieving a long service life for the battery pack!

In order to guarantee electrical functionality at all times, the bed should be connected to the mains as much as possible.

#### Safety cut-out device

The bed is equipped with an electrical, self-resetting safety cut-out device that prevents overloading of the motor systems. In the event of very severe overloading, the bed is automatically switched off.


# **Activation | Preparation**

#### Conditions for set-up

The bed is only approved for use in dry rooms (technical data sheet). A mains power supply is required for operation of the bed in any suitable room.

#### Hand control connection

The hand control should be connected to the socket provided, where appropriate. The spiral cable must be routed so that it is not under tension.



Routing the hand control cable

### **Bed transport**

The bed can be moved without auxiliary transportation devices. If necessary, release the brakes  $\square$  36).



**CAUTION** The bed should only be moved on a solid floor. Never attempt to push it over obstacles of over 2 cm in height. The maximum angle of inclination of the floor must not exceed 10°

**NOTE** The bed must be transported by at least two people, taking hold of the bed at the head and foot section.



Bed transport by at least two people

# **Activation | Electrical activation**



**WARNING** Ensure that the plug connection between the mains cable and the lying surface adjustment motor is correctly positioned!




**WARNING** Check the mains cable for damage. If the cable is damaged, the bed must not be used and must be taken out of service immediately!

**NOTE** Please be aware that incorrect handling of the mains isolation switch can lead to the battery pack (option) no longer charging.

Incorrect handling can include allowing the mains isolation switch to fall, pulling the mains isolation switch cable to release it from the socket and driving over the cable when transporting the bed.

1. Connect the mains isolation switch to the mains power socket.



2. Press the green button on the mains isolation switch for a second to allow the mains connection.

 Unlock the locking switch on the reverse of the hand control and the nurse keypad (option) using the key (open lock visible) in order to activate the bed's electrical functions.



# Activation | Using the battery pack (option) 1/2

The battery pack (option) allows the bed to be operated independently of the mains supply for at least two adjustment cycles.

The LED displays three colours:

Green	The battery pack is connected to the mains supply. The charging cycle is running.
Orange	The battery charger is being charged. The bed should not be operated independently of the mains supply.
Red	DANGER ZONE. The battery pack must be charged. The bed cannot be operated independently of the mains supply.
All lights off	The battery pack is fully charged: the mains cut-off is activated. No current flows in standby mode.

If you hear a beep, the battery pack needs to be recharged. The beep becomes weaker as the battery pack's charge diminishes. The battery pack is switched off shortly before deep discharging. After the bed is connected to the mains supply, press any button on the hand control to render it fully functional again. The battery pack is charged when it is connected to the mains after every use or if the charge has fallen too low.



**NOTE** If the bed is stored for a long period without being connected to the mains supply, the battery pack can discharge. The degree of discharge depends on the environmental conditions.

**NOTE** During the charging cycles, the battery pack is connected to the mains supply and is therefore supplied with electricity. The LED displays the battery pack's charging status during the charging cycle. The current cut-off is deactivated and current flows to the bed.



**WARNING** If electromagnetic interference occurs with other equipment in the area around the bed, please refrain from using these devices.

When the bed is being transported, it must always be handled carefully and protected from moisture.

V1.0 (06.07.2009)

# Activation | Using the battery pack (option) 2/2 and taking out of service

**NOTE** The bed is designed for use in an ambient temperature range from 10 °C to 40 °C, with a relative humidity of 30% to 40% and an air pressure of 700 to 1060 hPa.



**WARNING** The battery pack must only be replaced by personnel trained by Völker AG.



WARNING If the battery pack is faulty, degassing can occur. In rare cases, this can cause deformation of the battery pack housing. If this occurs, the bed must be immediately taken out of service and taken to an adequately ventilated room without sparks (electrical or fire). Immediately inform customer services should this occur!



WARNING The battery pack must be disposed of in an environmentally responsible manner using the appropriate services. Alternatively, you can return it to Völker AG for disposal.

To activate the hand control and the nurse keypad once the bed has been activated, the function key block must be disabled  $\square$  35.

### Taking out of service

The bed is taken out of service by disconnecting it from the mains supply. To do this, the mains plug is disconnected from the mains socket and, if necessary, the optional battery pack disconnected from the control unit. If the bed is to be out of service for a period longer than two weeks, then the 9V battery must also be disconnected from the motor.

To disconnect the 9V block battery, remove the two screws from the battery compartment lid on the two-side drive.



Then pull the contact strip off the battery.

# **Activation | Functional check**

### Visual inspection

Before each new occupancy of the bed, the following checks must be carried out:

- 1. Ensure the bed exhibits no visible signs of damage.
- Ensure that the insulation of the electrical cables is intact.
- Ensure that the next testing date has not been missed (see testing label).



**WARNING** Only undamaged beds that are still within their testing interval periods may be used!


#### **Functional check**

A functional check must be carried out before each new occupancy:

- All electrical functions must be actuated once to their terminal positions.
- The function of all side rails must be checked.
- 3. The braking function of the bed must be checked.

Once a fault-free functional check has been carried out, the bed is ready for use.

# Operation



# Operation | Key lock

Activating the key lock disables all of the bed's electrical functions.

The locking switches are located on the reverse of the hand control and the nurse keypad (option). They are unlocked and locked with the key (open lock visible).

If the bed functions cannot be actuated, check whether the key lock is activated.

The key should be removed from the bed when it is not required.



Hand control or keypad locked



Hand control or keypad unlocked



If the nurse keypad is locked, the hand control is also automatically blocked (system lock).

If only the hand control is locked, the nurse keypad retains its full functionality.

# **Operation | Central castor adjustment**

To brake the bed, step on the **red spot** on the double foot pedal. As soon as the double foot pedal is engaged in a 30° position, the bed is braked.

To move the bed into the desired position, move the double foot pedal into a horizontal position.

To angle the steering castors in the direction of travel, step on the **green spot** on the double foot pedal. As soon as the double foot pedal is engaged in a 30° position, a steering castor is fixed and the bed can be controlled safely.



#### **CAUTION "Risk of accident"**

If the bed is not being transported, the castors must always be in the braked position, since the bed may be required as a support for when the occupant/patient stands up or lies down. If the bed rolls away while unbraked, this can lead to a serious fall. After activating the central holding brake, check whether the bed is actually fixed, .i.e the castors are adequately braked.

The bed can be in a non-fullybraked status even after first activation or reactivation, and consequently it must be checked that the castors are correctly braked.



Bed braked



Everything free



Moving position (castors engaged)

# Operation | Side rails | General safety notes



**WARNING** All people whose duties involve manoeuvring of the side rails must have read and understood the following information:

- During actuation of the back, upper leg or lower leg section adjustment, of the lift or of the side rails, it must always be ensured that the occupant/patient is not touching the side rails, and that no part of their body is sticking out through the side rails.
- If the side rails are used for an individual whose psychological condition necessitates their use, it must be ensured that the hand control is stored completely out of their reach or its functions are completely locked. In all cases, care must be taken to ensure that entrapment risks are minimized
- Protective covers (cot side pads) are available as an accessory for the side rails. These provide additional protection against injury from contact with the side rails. The use of these protective covers is recommended for all persons for whom the risk of injury from unavoidable contact with the side rails is very high. Even with the covers, the care staff or occupant/patient must still take the necessary care when operating the bed.
- If the side rails are used, they
  must always be either completely raised and securely
  engaged, or completely lowered to the end stop. Because
  of the risk of entrapment, they
  must never be left in a position
  where they are not completely
  engaged.


 If the side rails are damaged, there is a risk that the patient will fall out of bed.

# Operation | Side rails | Model S 280

#### Raising the side rails

To raise the side rails, take hold of the grip recess and lift it first at the head end and then at the foot end until the rail audibly engages.





#### Lowering the side rails

Starting at the foot end, take hold of the recess grip and lift the top rail slightly. Press the button on the side and allow the side rail to lower. Repeat at the head end.



# WARNING "Risk of entrapment"

On model S 280, the side rails must be either completely raised and securely engaged or be in the fully lowered position.

For safety reasons, the side rails cannot be released when there is pressure being exerted on them from above.

Two configurations are possible for the side rails on model S 280:

All side rails down: the bed is accessible without restriction from both sides.

All side rails up: The occupant/patient enjoys maximum protection from rolling out of bed.

Contents 38 of 88 Help

# Operation | Side rails | Model S 380 1/2

#### Raising the side rails

- Pull the side rail element out horizontally until the end stop, and fold it upwards.
- 2. To adjust the height of the side rail, pull the telescopic section upwards until it reaches its end stop.

#### Lowering the side rails

 Press both buttons on the outside of the frame, right under the cross-member, to bring the heightadjustable rail element to its lowest position.



 Press the "Drücken / Push" trigger at the lower edge of the side rail element and tilt it sideways into the horizontal plane, so that it lies parallel to the floor.



3. Push the rail element completely underneath the lying surface.





The side rail elements can be used various times as required to protect the occupant of the bed/patient.

Raising all four rail elements offers the occupant/patient maximum protection.



WARNING The side rails should always be gripped with two hands at the ends of the element in question and guided upwards or downwards.



**CAUTION** Any weight on side rail elements that are pulled out horizontally must not exceed 15 kg!

# Operation | Side rails | Model S 380 2/2

**NOTE** Due to their exceptional stability, the side rail elements can also be used as a surface for storing bed linen (max. 15 kg) or as an additional supporting surface for care-related positions, such as Bobath treatments, or for the delivery of physiotherapy-related treatments.

Contents 40 of 88 Help

# Operation | Back section

The back lying surface can be adjusted using the hand control or the nurse keypad (option).

If necessary, disable the keypad lock 4 35.

The back section of the lying surface can be raised up to an angle of max. 70°.



**WARNING** When the back section is being raised with the side rails up, it must be ensured that none of the occupant's/patient's or any other person's body parts are sticking out through the side rails or are on top of them!



Back section up/down



Hand control



**WARNING "Risk of entrapment"** When adjusting the position of the back section, do not touch the frame in the area of the back section!



Back section up/down



Nurse keypad (option)



Contents 41 of 88 Help

# Operation | Mechanical rapid lowering of the back section / CPR function (option)

The S 280 and S 380 models can optionally be fitted with a mechanical rapid lowering function of the back section for resuscitation.



#### **WARNING**

The CPR function (**C**ardio**p**ulmonary **R**esuscitation function) may only be operated by qualified personnel!

Hold the back section on the mattress holder and pull the red lever located at the left or right underneath the seat section upwards to lower the back section rapidly. The back section can now be moved quickly downwards. The lowering process can be interrupted by letting go of the red lever.



**WARNING** The back section must always be held on the mattress holder so that the bed is not lowered suddenly with the occupant/ patient.

**NOTE** When the CPR function has been used, the back section must then always be lowered to its lowest position. This can be done by activating the "back

section down" button on the hand control or the nurse keypad or by activating the red lever on the left or right underneath the seat section again (CPR function) and lowering the back section to its lowest position.

If this note is not observed, the back section adjustment range of 70° can no longer be fully utilised.

Red lever for mechanical rapid lowering of the back section for resuscitation:



# **Operation | Upper leg section**

The position of the upper leg lying surface can be adjusted using the hand control or the nurse keypad (option).

If necessary, disable the keypad lock  $\square$  35.

The upper leg section of the lying surface can be raised up to an angle of max. 45°.



**WARNING** When the upper leg section is being raised with the side rails up, it must be ensured that none of the occupant's/ patient's or any other person's body parts are sticking out through the side rails or are on top of them!

Please note that the position of the lower leg section can be adjusted by pulling on the mattress holder (or on model S 280 on the lying surface) (see also next page).



**WARNING** "Risk of entrapment with the S 380" When the position of the upper leg section is being adjusted, there is a risk of entrapment between the raised side rail and the foot section.



Upper leg section up / down



Hand control



Upper leg section up / down



Nurse keypad (option)

# Operation | Lower leg section

The lower leg section can be moved manually to any position of maximum 45° angle by pulling on the mattress holder (or the lying surface on model S 280).

To lower the lower leg section, the mattress holder (or the lying surface on model S 280) is pulled and the section lifted to the end stop and then lowered. The locking mechanism is disengaged automatically.



**WARNING** When the lower leg section is being raised with the side rails up, it must be ensured that none of the occupant's/patient's or any other person's body parts are sticking out through the side rails or are on top of them!





# **Operation | Lying surface height**

The position of the entire lying surface can be adjusted using the hand control or the nurse keypad (option).

If necessary, disable the keypad lock 🕮 35.

The lying surface height can be adjusted from 40 cm to 80 cm.



**WARNING** To avoid dangers to the occupant/patient from falling, we recommend that the bed be lowered all the way except when delivering care.



**DANGER** Before lowering the bed, it must be ensured that no people, limbs, pets, bed linen or other objects are trapped between the lying surface and the lower frame or floor. The bed's position must be stable (castors braked) when the patient is getting in and out of the bed!



**WARNING** When the height adjustment mechanism is actuated with the side rails up, it must be ensured that the occupant/ patient does not have any contact with the side rails, or that any parts of either the patient's body or of other persons are sticking through the side rails!



Lying surface up / down



Hand control



Lying surface up / down



Nurse keypad (option)

# Operation | Anti-Trendelenburg<sup>1</sup> and Trendelenburg positioning<sup>2</sup>

Trendelenburg positioning can only be set using the nurse keypad (option).

If necessary, disable the keypad lock 🕮 35.

The Trendelenburg position can be adjusted up to an angle of 12°.



**CAUTION** If a fault occurs with the lifting function, or the mains power supply fails and the battery pack is completely discharged, the Trendelenburg function cannot be engaged. The occupant/patient may have to be placed in a different bed. With the optionally-available battery pack, all of the bed's functions remain available to you even in the event of a power failure.



**WARNING** Since the Trendelenburg position depends on clinical indications, it must only be used with appropriate clinical approval.

The lying surface is automatically returned to its horizontal position if it is moved to its highest or lowest position.



Anti-Trendelenburg positioning <sup>1</sup>



Hand control



Anti-Trendelenburg positioning <sup>1</sup> and Trendelenburg positioning <sup>2</sup>



Nurse keypad (option)

<sup>&</sup>lt;sup>1</sup> Head end raised

<sup>&</sup>lt;sup>2</sup> Head end lowered

# **Operation | Comfort seating position**



Occupants/patients who are unable to leave the bed either because their circulatory system is too unstable or they first have to be "taught" how to sit can gain a lot of benefit from the comfort seating position as it provides an active sitting position in the bed.

Anti-Trendelenburg positioning



Auto-Contour

Hand control with anti-Trendelenburg positioning and Auto-Contour (option)

#### Adjusting the comfort seating position

1. Move the back



and the upper leg section upwards a little to a comfortable position.

You can alternatively reach this position in a single step by pressing

the Auto-Contour button



2. Swivel the bed to the comfort seating position by pressing the Anti-Trendelenburg button.

#### Restoring the straight lying surface

To return to a horizontal lying position, move the lying surface -- and the back and upper leg section in any preferred order to their lowest position.

**Contents** 47 of 88 Help

# **Operation | Bed extension (option)**

#### Model S 380

The bed can be extended using a telescopic bed extension (option) by 20 cm.

To pull out the bed extension, pull the two pins located on the underside of the bed extension simultaneously downwards.





# Cleaning and disinfection



Contents 49 of 88 Help

# Cleaning and disinfection 1/2

In order to maintain consistent functioning, the care bed should be cleaned, disinfected and tested at the soonest possible time following each use, so that it can be reused immediately without risk. Incorrect cleaning/disinfection of the bed can cause danger.

#### Cleaning

Depending on the degree of soiling, we recommend cleaning the bed with a damp cloth or similar.

#### Wipe and spray disinfection

For wipe and spray disinfection, the disinfecting agents featured in the VAH (Verbund für angewandte Hygiene e.V.) list dated 01.01.2008 can be used in their specified concentrations. They must be applied at the dilution ratio specified by the relevant manufacturer's instructions for use. The list can be obtained under ISBN Number 978-3-88681-089-5.

**NOTE** Solvents are not permitted. Scouring agents, abrasive sponges or other blunting agents must not be used.

Organic solvents such as halogenated / aromatised hydrocarbons and ketones must not be used.

Depending on the degree of soiling, we recommend cleaning the bed with a damp cloth or similar.

The following instructions must be noted when using cleaning and disinfecting agents:

- The working solutions should generally be freshly prepared.
- The concentrations used should be neither higher nor lower than those given in the list. The so-called "shot" method should not be used under any circumstances. Under no circumstances should someone using a disinfectant follow their own judgement to add a detergent such as a soap or detergent substance (leads to soap failure).

- There is a risk of fire or explosion from alcoholic spray disinfectants when these are used over large areas.
- They must not contain any corrosive or irritant components.
- They must not contain any substances that change the surface structure or gripping properties of the materials.
- · Lubricants must not be touched.
- The pH value of the water must be no higher or lower than 6 - 8.
- Water should not exceed a total water hardness of 0.9 mmol/l (up to 5 deg d). (Desalinated water should not be used).

# Cleaning and disinfection 2/2

 $\begin{array}{lll} \text{Chloride} & < 100 \text{ ppm} \\ \text{Silicates as SiO}_2 & < 15 \text{ ppm} \\ \text{Iron} & < 0.05 \text{ ppm} \\ \text{Manganese} & < 0.01 \text{ ppm} \\ \text{Copper} & < 0.05 \text{ ppm} \end{array}$ 

The specifications in the VAH list, as well as the specifications we have issued do not absolve the user of the obligation to carry out his or her own checks and investigations, since the ratios (e.g. water hardness) can vary depending on the geographical location. It is therefore impossible to offer legally binding assurance of certain properties.

If unsuitable washing and disinfecting agents are used, if the mixing ratio is incorrect or if there is inadequate care of the beds, damage can occur to the surface coating for which Völker AG is not liable.



# WARNING "Risk of electric shock / fire and functional failure"

The bed must always be disconnected from the power supply during cleaning and disinfection.

The plug and the socket of the hand control are only protected against spray when the control is stowed and the appropriate cover is in place.

# Spray lances and automatic washing systems

Cleaning and disinfection using spray lances from high-pressure cleaning equipment and in automatic bed washing systems is **not** permissible.

#### **Maintenance**

The **Maintenance** section contains information on how to carry out maintenance work.

#### **CONTENTS**

Staff training 4 53

Safety notes 4 53

Maintenance schedule 4 54



# Maintenance | Staff training, safety notes

#### Staff training

Every person involved with maintenance or servicing must at least have read

- · the safety notes and
- the service manual

and be qualified in accordance with MPBetreibV §§ 4 and 6. To avoid errors and ensure that the care beds work properly, these documents must always be accessible to service personnel.

Before commencing maintenance work, the service manual and the instructions for use must be read in full by the person responsible for carrying out the servicing.

#### Safety notes

During maintenance and technical checks, the following specifications must be strictly observed:

- The room's electrical installations must satisfy the requirements of the current state of the technological art and the care bed must be used accordingly.
- The castors must be placed in the "brake position".
- The care beds are not protected against explosion and must therefore only be maintained in an environment free from flammable substances and materials.



**WARNING** Before carrying out any maintenance work, always disconnect the mains plug **and** disconnect the bed from the battery pack.



**WARNING** Maintenance and repair work should only be carried out after the bed has been disinfected



WARNING After maintenance (repair) work is complete, always check that the bed is functioning correctly. It must be checked that the bed can be used correctly without risk to the occupant/patient or care personnel.

# Maintenance | Maintenance schedule

The care beds require little maintenance. All movable parts for the height adjustment mechanism, the lying surface motor systems and the side rails are provided with long-lasting lubricant at the factory. It is recommended that the care beds are subjected to a regular, or at least once a year, technical check 55 (incl. visual inspection and functional check) as described in the checklist \$\mathbb{\text{\text{\text{u}}}}\$ 80 and any damage uncovered as a result, such as signs of wear and tear, loose screws or breaks/fractures, be eliminated immediately.

After every lengthy period of non-use, a technical check 55 must be carried out.

Period	Work to be carried out
Annually and after lengthy periods of non-use	Technical check 🕮 55
Every 2 years*	Replace the 9V block battery (see below)
As required	Lubricate mechanical parts Replace the battery pack or other wear parts (see service manual)

#### Replace the battery

To replace the 9V block battery, remove the two screws from the battery housing on the double drive.


Use only brand-name batteries, and dispose of spent batteries in an environmentally friendly way.



\* as well as after every emergency lowering of the back section with a 9V block battery, if the bed does not have a battery pack (option).

#### Technical check

The section entitled **Technical check** contains all the information needed to carry out the technical check in accordance with MPBetreibV, BGVA3, UVV on care beds and measurement as per VDE 0751-1. Other (e.g. country-specific) specifications have not been included here. This does not absolve the operator from the obligation to observe any such specifications.



#### **CONTENTS**

Checking the trapeze bar grab handle

Other accessories 

58

**58** 

#### Technical check 1/3

#### 1. Visual inspection


Check the frame parts for plastic deformations and / or wear and tear. These include the lower frame, the lifting mechanism, all parts of the lying surface (back, seat, upper leg and lower leg sections), trapeze bar, trapeze bar holder and castors.


#### 2. Functional check of the side rails

Check that the locking mechanism for the side rails is working correctly and that there are no visible deformations or signs of wear and tear on the side rails.

Check that the prescribed distances are maintained, even when the side rails are placed under load. All measurements of side rail distances must be carried out in the flat lying surface position.

The measurement of A must be performed using a cone as defined in IEC 60601-2-38. The force used for measurement A must be at least 50 N.





Letter	Description	Dimension
Α	The greatest distance between the elements within the scope of the SIDE RAIL in a raised / engaged position or the area formed by the SIDE RAILS and fixed parts of the BED.	≤120 mm
В	Thickness of the mattress for CORRECT USE	defined by the manufacturer
С	Height of the upper edge of the SIDE RAILS above the mattress (see "B") without compression	≥220 mm
D	Distance between the HEAD or FOOT SECTION and the SIDE RAIL	≤60 mm or ≥235 mm
E	Distance between the separated SIDE RAIL with the LYING SURFACE in flat position	≤60 mm or ≥235 mm
F	Distance of all accessible openings between the SIDE RAILS and the LYING SURFACE	if D or E E≥235 mm ≥60 mm then F ≤60 then F mm ≤120 mm
G	Total length of the SIDE RAIL or sum of the length of the divided SIDE RAILS on one side of the BED	∑ Gx ≥1/2 the length of the LYING SURFACE

Contents 56 of 88 Help

#### Technical check 2/3

#### 3. Functional check of the brakes

Check that the brake is functioning correctly (braked, steer, free running).

#### 4. Functional check of the motors

Travel through the full adjustment range of each motor. Look out for any unusual noises, watching the speed, ease of running etc., and check that the selected function travels in the correct direction. Particularly ensure that the motor switches off automatically when it reaches its terminal position\* (terminal position switch).

\* Please note that the terminal position can vary from bed model to bed model. Please consult the technical specifications for this, or if in doubt, contact our service department.

#### 5. Mains connection cable

Check

- the mains connection cable, incl. cable guides,
- the strain relief, including kink protection sleeve,
- the mains connection plug for damage.

#### 6. Cabling

Check the cable guides and that the plug connectors are correctly seated and do not exhibit any damage.

Check the cables for damage.

#### 7. Housing

Check all housings for damage. All screws must be firmly tightened and seals must not exhibit any visible damage.

#### 8. Mechanical check

Check the function of the Rastomat or the pneumatic spring (where fitted) by manually moving the lower leg section to the individual positions.

#### 9. Changing the battery

Replace the 9V block battery at 2-yearly intervals. Only use brand-name batteries and dispose of the old ones in an environmentally responsible manner (Replacing the batteries 4).

#### 10. Measurement as per VDE 0751-1

The bed must be checked electrically in accordance with VDE 0751-1. The equivalent unit current leakage must be measured.

### Technical check 3/3

#### 11. Trapeze bar grab handle


Check whether the plastic and holding frame of the grab handle exhibit any damage and that the fixing rods on the trapeze bar are intact.

The handle of the trapeze bar, including the fastening strap, must be replaced every five years.

#### 12. Other accessories

Other accessories must always be checked in accordance with the manufacturer's instructions.

# **Troubleshooting**



Contents 59 of 88 Help

# Troubleshooting | Table of faults 1/2

The following table contains information on possible functional problems that the users can resolve themselves.

Faults that can only be eliminated by technical engineers are described in the service manual.

Anyone responsible for carrying out maintenance and repair work must at least have read the safety notes and the service manual and be qualified in accordance with MPBetreibV §§ 4 and 6.

**NOTE** Before carrying out any troubleshooting, ensure that the battery pack (option) is charged (the yellow LED flashes during the charging process at intervals commensurate with the charge level) and that the bed is connected to the mains supply (the mains plug is in a live socket).




warning Before carrying out any repair work, ensure that the bed has been disconnected from the power supply and the battery pack has been disconnected if present.

# Troubleshooting | Table of faults 2/2

Fault	Possible cause	Troubleshooting
The bed cannot be adjusted electrically.	The keypad lock is engaged.	(1) Deactivate the key lock 🚨 35.
		(2) Connect the plug or check the socket.
	<ul><li>socket is not live.</li><li>The battery pack (options) is not con-</li></ul>	(3) Press the green button on the mains isolation switch¹ while simultaneously actuating
	nected or is flat.	any function on the hand control.
	The hand control is faulty.	(4) Replace the hand control.

<sup>1</sup> Green button on the mains isolation switch:



The mains isolation switch ensures that when no electrical function is being actuated, there is no mains voltage to the bed. (Exception: while the battery pack (option) is being charged, there is mains voltage to the bed. This is indicated by a flashing LED on the battery pack).

# **Troubleshooting | Service points**

If necessary, please seek assistance from the relevant contact at your nearest Völker sales organisation. You will receive all the information you need for comprehensive service promptly.

Contents 62 of 88 Help

# **Appendix**

The **Appendix** section contains the technical specifications and classifications, details on the service life and disposal of the equipment and links to the manufacturers' declarations and forms found in the appendix.

Konformitätserklärung Anhang VII EU-Richtlinie 93/42/EWG

Der Unterzeichnende Völker AG Wullener Feld 79 58454 Witten

mit einer Fertigungsstätte unter der Adresse Völker AG Ahornstraße 4

bestätigt, dass die nachfolgend bezeichneten Produkte in der von uns in Verkehr gebrachten Ausführung die grundlegenden Anforderungen des Anhangs I der EU-Richtlinie 93/42/EWG erfüllen. Es wurden die folgenden Normen angewendet:

\_ ... \_ . . . . . . . . .

Declaration of conformity
Appendix VII
EU Directive 93/42/EEC

The signatory Völker AG Wullener Feld 79 58454 Witten/Germany

with a manufacturing side at: Völker AG Ahornstraße 4 09661 Hainichen/Germany

confirms that the products described below and in the form distributed by purselves meet the basic requirements of Appendix I of EU Directive 93/42/EEC. The following standards are applied:

DIN EN 60601-1,

Déclaration de conformité Annexe VII Directive EU 93/42/CEE

La soussignée Völker AG Wullener Feld 79 58454 Witten/Allemagne

avec un site de production à: Völker AG Ahornstraße 4 19661 Hainichen/Allemagne

confirme que les produits spécifies cidessous sont conformes, dans le modèle mis en circulation, aux exigences fondamentales de L'annexe L de la directive européenne 93/42/CEE. Les standards suivants sont appliqués:

----

#### **CONTENTS**

Symbols used 🕮 64

Technical data 🕮 65

Type label 🕮 68

Service life / disposal 🕮 70

Manufacturer's declarations, forms, electronic instructions for use 4 71

Contents 63 of 88 Help

# Appendix | Symbols used



#### Warning symbols

Information marked with this symbol must always be read and strictly observed!



Warning about the risk of crushing and entrapment!



Note the information in the instructions for use!



Safe working load 230 kg



Max. permissible weight of the occupant/patient 200 kg (approx. 30 kg for the mattress, trapeze bar and accessories)



Direct current



Alternating current



Protection class II device, protection-insulated



Type B device as per DIN EN 60601-1



The product satisfies the fundamental requirements of Appendix 1 of EU Directive 93/42/EEC.



The product must be disposed of in accordance with

EU Directive 2002/96 pertaining to old electrical and electronic equipment.

<sup>\*</sup> this value is 20 kg lower for beds built before 10/2008

# Appendix | Technical data (standard design) 1/2

The details marked below with \* are dependent on the length, width or model of the bed. The values specified relate to model 380 S.

External dimensions (W x L)	
for lying surface 90 x 200 cm *	98 x 210 cm
for lying surface 90 x 210 cm *	98 x 219 cm
for lying surface 90 x 220 cm *	98 x 229 cm
for lying surface 100 x 200 cm *	109 x 210 cm
for lying surface 100 x 210 cm *	109 x 219 cm
for lying surface 100 x 220 cm *	109 x 229 cm
for lying surface 110 x 200 cm *	119 x 210 cm
for lying surface 110 x 210 cm *	119 x 219 cm
for lying surface 110 x 220 cm *	119 x 229 cm
Height of the lower edge (min./max.) *	approx. 25 cm / 65 cm
Height of the upper edge (min./max.) *	approx. 76 cm / 116 cm
Height of the upper edge of the lying surface *	approx. 40 cm - 80 cm
Lying surface (4-part) *	90 x 200/210/220 cm
	100 x 200/210/220 cm
	110 x 200/210/220 cm
Volumetric weight of the	40 - 50 kg/m <sup>3</sup>
mattress material	-
Dead weight *	135 kg

Safe working load for the bed	230 kg** As a result of the weight (approx. 30 kg) of the mattress, trapeze bar and other accessories, the maximum permissible weight of the occupant/ patient is 200 kg. If another mattress is used or other accessories, this value must be recalcu- lated!
Safe working load for the trapeze bar holder	75 kg
Safe working load for drip stands	2 kg / hook
Castors	4 pcs, Ø 150 mm
Max. castor load	100 kg (dynamic)
Mains voltage Rated power Rated frequency	AC 220 - 240 V, 100 - 120 V 173 VA 50 / 60 Hz

Contents 65 of 88 Help

<sup>\*\*</sup> this value is 20 kg lower for beds built before 10/2008

# Appendix | Technical data (standard design) 2/2


Primary fuse	2.0 A
Hand control fuse	Type: Polyswitch RXE 025
Lying surface motor fuse	Type: Polyswitch, solid, 2.5 A
Lifting motor fuse	Type: Polyswitch, solid, 3.75 A per motor
Battery	Type: 9 V block battery (alkali manganese, commercially available)
Battery pack (option)	Type: 4 x 6 V block battery (lead gel) 1.2 Ah
Operating temperature range	+ 10 °C to + 40 °C
Transport / storage tempera- ture range	- 20 °C to + 60 °C
Air humidity	30 % to 75 % rel.
Atmosphere range	700 hPa to 1060 hPa

# Appendix | Classification

Protection against electric shock	Protection class II or device with internal electrical current sour
Protection type by housing as per EN 60259	IPX4 not suitable for cleaning in automatic washing systems
Degree of protection of the applied part against electric shock as per DIN EN 60601-1	Type B
Degree of protection against explosive substances and mixtures	The bed is not protected against explosion and should not be used in an environment in which flammable anaesthetics or flammable cleaning agents are present (see brochure from accident prevention organisation ZH 1/200)
MPG grouping	Class I

Mode of operation	Int. 1 min. / 10 min. On-time max. 1 min. Off-time 10 min.
Technical check	1x yearly

# Appendix | Type label 1/2




The type label is located on the inside of the head panel.

Raise the back section to read the type label.



Type specifications	Declaration
1st line	Model designation. In the example: NURSING BED Model S 380
2nd line	ID No. (made up as follows): S 380 = Model -2009 = Year of construction -7 = Production week
	AUS = Mains plug version (e.g. AUS = Australia)
3rd line	Input: mains voltage; mains frequency; consumed power
4th line	Physical life time: Max. uninterrupted on-time of electric motor adjustment. In the example: 1 min./10 min. In other words, the bed may be operated with the electric motors for max. 1 min uninterrupted within a 10-minute period (protection against overheating).
5th line	Safe working load (accessories, mattress and weight of occupant/patient).
6th line	Suitability for automatic washing systems. In the example: Unsuitable with automatic washing system.

# Appendix | Type label 2/2



Measured values Declaration

1st line Equivalent unit current leakage in µA

The specified initial measured values were measured in accordance with VDE 0751-1.

The barcode (code 39) contains the numeric ID No. (10 digits).

# Appendix | Service life / disposal

The care bed's expected service life is approx. 10 years. To ensure environmentally responsible disposal after decommissioning, please contact your responsible area sales manager.

Contents 70 of 88 Help

## Appendix | Manufacturer's declarations, forms, electronic instructions for use

#### Manufacturer's declarations

- Guidelines and manufacturer's declaration Electromagnetic compatibility (6.8.3.201 a) 6)) 74
- Table 204 Guidelines and manufacturer's declaration Electromagnetic compatibility for all devices and systems that are not life-sustaining (6.8.3.201 b)) 76
- Table 206 Recommended protected distance between portable and mobile HF telecommunications equipment and the S 280 or S 380 – for devices and systems that are <u>not</u> life-sustaining (6.8.3.201 b)) 78

#### **Forms**

- Technical check as per MPBetreibV, BGVA 3, UVV on hospital and care beds, incl. measurement as per VDE 0751-1 80

#### Electronic instructions for use

- Requirement for the use of the electronic instructions for use 4 82
- CD-ROM with electronic instructions for use \$\omega\$ 83

Contents 71 of 88 Help



# Konformitätserklärung / Declaration of Conformity / Déclaration de conformité

Konformitätserklärung Anhang VII EU-Richtlinie 93/42/EWG

EU-Richtlinie 93/42/EWG Der Unterzeichnende Völker AG

Wullener Feld 79

58454 Witten

mit einer Fertigungsstätte unter der Adresse: Völker AG Ahornstraße 4 09661 Hainichen bestätigt, dass die nachfolgend bezeichneten Produkte in der von uns in Verkehr gebrachten Ausführung die grundlegenden Anfordeungen des Anhangs I der EU-Richtline 93/42/EWG erfüllen. Es wurden die fölgenden Normen angewendet:

DIN EN 60601-1, DIN EN 60601-1-2, DIN EN 60601-2-38 (Anwendungen der relevanten Teile), EN 1970 (Anwendungen der relevanten Damit sind die Anforderungen des Medizinproduktegesetzes zur Anbingung einer CE Kennzeichnung erfüllt.

Bei einer nicht mit dem Hersteller abgestimmten Änderung des Produktes verliert diese Konformtiätserklärung ihre Gütigkert.

Bezeichnung der Produkte : Pflegebetten 5380, 2080, 3080, S 380 und S 280. EG-Richtlinien:
Richtlinie 93/42/EWG vom 14.06.1993 über
Medizinproudkte (Anhang I., Grundlegende Anfordeungen?).
Die Produkte sind Produkte der Klasse I gemäß Anhang VII des Medizinproduktegesetzes MPG vom Medizinproduktegesetzes MPG vom

To print out this document, please use the

Witten, 22.01.2009

CD-ROM.

corresponding

PDF file on the

Declaration of conformity

Appendix VII EU Directive 93/42/EEC

The signatory Volker AG Willener Feld 79 58454 Witten/Cermany

with a manufacturing side at: Völker AG Ahornstraße 4 09661 Hainichen/Germany confirms that the products described below and in the form distributed by ourselves meet the basic requirements of Appendix for EU Directive 93/42/EEC. The following standards are applied:

DIN EN 60601-1,
DIN EN 60601-1-2,
DIN EN 60601-2-38 (applicable parts
only),

EN 1970 (applicable parts only).

DIN EN 60601-1

The requirements of the medical products law pertaining to the display of a **CE seal** of approval are thereby fulfilled.

This declaration of conformity becomes invalid if the products are altered without the agreement of the manufacturer.

Description of products Type/Article No.: Nursing beds 5380, 2080, 3080, S 380 and S 280. EU Directives:
Directive 33/42/EEC of 14.06.1993
concering medical products (Appendix I, Basic requirements). The design and construction of this product conforms to Class I (Appendix VII) Medical products law (MPG) of 02.08.1994.

Déclaration de conformité Annexe VII Directive EU 93/42/CEE

La soussignée Völker AG Wullener Feld 79 58454 Witten/Allemagne avec un site de production à: Völker AG Ahornstraße 4 09661 Hainichen/Allemagne confirme que les produits spécifies cidessous sont conformes, dans le modèle mis en circulation, aux exigences iondamentales de L'annexe L de la afriective européenne 93/42/0E. Les standards suivants sont appliqués :

DIN EN 60601-1-2, DIN EN 60601-2-38 (en partie applicable), EN 1970 (en partie applicable). Les exigences de la loi sur les pn

Les exigences de la loi sur les produits médicaux concemant leport de la marque CE sont ainsi satisfaites.

marque CE sont ainsi satisfaites. Cette déclaration de conformité est Désignation des produits Modèle/Référence : Lits médicalisés 5380, 2080, 3080,

S 380 et S 280.

produits, non autorisée par le fabricant.

nvaledée en cas de modification des

Directives européennes :
Directive 9342/CEE du 14.06.1993 sur
les produits médicaux (annexe l
« Exigences fondamentales »). La conception du produit est conforme à la classe I (annexe VIII). Loi sur les produits médicaux (MPG) du 02.08.1994.

Info

Heinrich Volker Vorstandsvorsitzender / Executive board (chair) / Dirctoire (Président)

# Table 201 - Guidelines and manufacturer's declaration - Electromagnetic compatibility (6.8.3.201 a) 3))

The S 280 and S 380 are intended for use in the electromagnetic environment as described below. The customer or user of the bed should ensure that it is used in such an environment.

Emissions measurements	Conformities	Electromagnetic environment - guideline				
HF emissions IEC 61000-3-2	Class A	The S 280 and S 380 are suitable for use in all institutions, including residential areas and such that are directly connected to a public mains supply network that also supplies buildings intended				
Fluctuations in voltage / flicker	Satisfied	for residential purposes.				
IEC 61000-3-3						
RF emissions	Satisfied	The S 280 and S 380 is not suitable for connection to other equip-				
CISPR 14 – 1		ment.				

### Guidelines and manufacturer's declaration - Electromagnetic compatibility (6.8.3.201 a) 6))

The S 280 and S 380 are intended for use in the electromagnetic environment as described below. The customer or user of the bed should ensure that it is used in such an environment.

Compatibility check	IEC 60601 Testing level	Conformity level	Electromagnetic environment - guideline
Discharge of static electricity (ESD) IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	The ground should be made of wood or concrete and be laid with ceramic tiles. If the ground is covered with synthetic material, the relative air humidity must be at least 30%.
Rapid transient electrical disturbance / bursts IEC 61000-4-4	± 2 kV for voltage supply ± 1 kV for input and out- put cables	± 2 kV for voltage supply Not suitable!	The quality of the voltage supply should be commensurate with a typical commercial or hospital environment.
Surge voltages IEC 61000-4-5	± 1 kV series mode voltage ± 2 kV common-mode voltage	± 1 kV series mode voltage  Not suitable!	The quality of the voltage supply should be commensurate with a typical commercial or hospital environment.

© Continued on next page.

Compatibility check	IEC 60601 Testing level	Conformity level	Electromagnetic environment - guideline			
Voltage drops, short interruptions and voltage fluctuations in the mains voltage supply	$< 5 \% \ U_{\scriptscriptstyle T}$ (>95 % voltage peaks in $U_{\scriptscriptstyle T}$ ) for 0.5 cycles	$< 5 \% U_{\rm T}$ (>95 % voltage peaks in $U_{\rm T}$ ) for 0.5 cycles	The quality of the voltage supply should be commensurate with a typical commercial or hospital environment.  If the user requires continued function			
IEC 61000-4-11	40 % U <sub>T</sub> (60 % dip in U <sub>T</sub> ) for 5 cycles	40 % U <sub>T</sub> (60 % dip in U <sub>T</sub> ) for 5 cycles	from the S 280 or S 380 when inter- ruptions in the power supply occur, it is recommended that the S 280 or S 380 be powered via an uninterruptible power			
	$70 \% U_{\scriptscriptstyle T}$ (30 % dip in $U_{\scriptscriptstyle T}$ ) for 25 cycles	70 % U <sub>T</sub> (30 % dip in UT ) for 25 cycles	supply or a battery.			
	$< 5 \% \ U_{\scriptscriptstyle T}$ (>95 % dip in $U_{\scriptscriptstyle T}$ ) for 5 sec	$< 5~\%~\mathrm{U_T}$ (>95 % dip in $\mathrm{U_T}$ ) for 5 sec				
Magnetic field at supply frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at mains frequency should correspond to the typical values present in commercial and hospital environments.			
Remark 1: U <sub>T</sub> is the mains AC voltage before the testing level is applied.						

## Table 204 – Guidelines and manufacturer's declaration – Electromagnetic compatibility for all devices and systems that are <u>not</u> life-sustaining (6.8.3.201 b))

The S 280 and S 380 are intended for use in the electromagnetic environment as described below. The customer or user of the bed should ensure that it is used in such an environment.

Compatibility check	IEC 60601 Testing level	Conformity level	Electromagnetic environment - guideline
			Portable and mobile radio equipment should not be used any closer to the S 280 or S 380, or their cables, than the recommended protective distance which is calculated based on the equation applicable for the transmission frequency.
			Recommended protective distance
Conducted HF	3 Vrms	3 V	$d = 1.17\sqrt{P}$
disturbances IEC 61000-4-6	150 kHz to 80 MHz		$d = 1.17\sqrt{P}$ 80 MHz to 800 MHz
			$d = 2.33\sqrt{P}$ 800 MHz to 2.5 GHz

© Continued on next page.

Compatibility check	IEC 60601 Testing level	Conformity level	Electromagnetic environment - guideline
Radiated HF disturbances	3 V/m 80 MHz to 2.5 GHz	3 V/m	Where <i>P</i> is the rated output of the transmitter in Watts (W) as specified by the transmitter manufacturer and <i>d</i> is the recommended protective distance in metres (m).
120 01000 10			The field strength of stationary radio transmitters should be less than the conformity level <sup>b</sup> at all frequencies, as verified by an on-site test <sup>a</sup> .
			Disturbance is possible in the environment of equipment that bears the following label.

Remark 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Remark 2: These guidelines may not be applicable in all cases. The spread of electromagnetic parameters is influenced by the absorption and reflection of buildings, objects and people.

<sup>a</sup> The field strength of stationary transmitters, such as base stations of cordless telephones and mobile land radio, amateur radio stations, AM and FM radio and TV stations cannot be quantified in theory in advance accurately. To determine the electromagnetic environment in relation to a stationary transmitter, a study of the site should be considered. If the measured field strength at the site where the S 280 or S 380 is used exceeds the conformity levels stated above, the S 280 or S 380 should be observed to ensure that it functions correctly. If unusual features are observed, additional measures may be required, such as modified alignment or a different location for the S 280 or S 380.

In the frequency range from 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

Contents 77 of 88 Help

Table 206 - Recommended protective distance between portable and mobile HF telecommunications equipment and the S 280 or S 380 – for devices and systems that are <u>not</u> life-sustaining (6.8.3.201 b))

The S 280 or S 380 are intended for operation in an electromagnetic environment in which HF disturbances are controlled. The customer or user of the S 280 or S 380 can help to avoid electromagnetic disturbances by respecting the minimum distance between portable and mobile HF telecommunications equipment (transmitters) and the S 280 or S 380 The recommended minimum distance d is dependent on the maximum power output of the communication device (see below).

Use in environments specified in which the radiated RF disturbances are controlled. The buyer or user of the S 280 or S 380 can help avoid electromagnetic interference by respecting the minimum distance between portable and mobile RF communications equipment and the S 280 or S 380. The recommended minimum distance is dependent on the maximum power output of the communication device.

	Protective distance according to the frequency of the transmitter m				
Rate output of the transmitter	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz		
W	d = 1.17 P	d = 1.17 P	d = 2.33 P		
0.01 0.12 √		0.12 √	0.23 1		
0.1	0.37	0.37	0.74		
1	1.17	1.17	2.33		
10	3.69	3.69	7.38		
100	11.67	11.67	23.33		

© Continued on next page.

For transmitters whose maximum rated output is not specified in the table above, the recommended protective distance d in metres (m) can be calculated using the equation in the column that features P as the maximum power output of the transmitter in Watts (W) as defined by the transmitter manufacturer.

Remark 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Remark 2: These guidelines may not be applicable in all cases. The spread of electromagnetic parameters is influenced by the absorption and reflection of buildings, objects and people.



Technical check of Völker hospital and healthcare beds in accordance to German standards and safety regulations

incl. measurements required

Proje	Project/customer no.:	Location of the bed:	Type o	Type of bed/product :	nct :		
Identi	Identification of the hospital and healthcare beds (e.g. due to the Inventory	beds (e.g. due to the inventory					
of the Date	of the facility and Volker bed identification Date of check:	number);	Name of technician:				
				ĺ			
Kind	Kind of check	Component to be checked		Anually	Accepted	Not accepted	Not applicable
Visua	Visual inspection	Inscription on device readable					
		Instructions for use available		ż			
		Lying surface		a å			
		Trapeze bar adapter/infusion bar adapter	adapter				
		Power supply cable and plug/charging connection (charging connection part of 980 S 980.1 and S 980.1 MIT)	ging connection	·B			
		Pull relieve/bend protection		å			
		Fit of plug-and-socket connections	9	B			
		Motor housing/transformer housing	6	B.			
		Housing of hand control		B			
		Cable of hand control		e B			
		Nurse control		å.			
		Trapeze bar, assist rall spacer, additional accessories	ditional accessories	B*Æ*			
		Stroke head traverse, head and foot section	oot section	å å			
		Castors  Desired of the continuers					
		Side rails, if applicable with telescope	ado	. å			
		HilLow-elevation: check screw locking (only for 5380)	king (only for 5380)	°S			
Funct	Functional inspection of side rails,	Locking devices		*			
if app.	licable with telescope	Deformation		×			
		Abrasions		*×			
Funct with h conso	Functional inspection of drives with hand control and caregiver's console	Lying surface, back section, upper leg section, lower leg section, hely adjustment, Trendeenburg position / reverse Trendeenburg position, base frame, approach all end switches	r leg section, lower leg elenburg position / reverse ches	x*/M*			
		Angle limitation of back section and upper leg section > 90° (only for hospital beds)	nd upper leg section > 90°	*×			
		Catch adjustment of foot section/gas spring	gas spring	×			
		CPR function (if existing)		×			
		Brake (electrical of mechanical)	- brake applied	*×			
			- brake not applied				
		(only for hospital beds and S 280 / S 310 / S 380)	- straight forward				
		Mechanical release (only for electrical brakes of hospital beds)	rical brakes of hospital	×			
Funct	Functional inspection replacement	9V battery (only for healthcare beds and S 960-2) Replaced (Yes/No)	ds and S 960-2)	. A2*			
		Trapeze bar handle and belt (if existing) Replaced (Yes/No)	isting)	A5*			
Funct	Functional inspection miscellaneous	Bed extension (if existing)		B.			
		Bedding storage / bedding drawer (if existing)	(if existing)	å			
Comment	nent						
JOL	Leakage current < 5000 µA, in accordance to VDE 0751			γи			
Sgally relevant	Potential equalization impedance < 0.2 Ohm, in accordance to VDE 0751 (if existing)			а			
п	Measuring instrument S/N						
Total	Total result of the inspection:						
Signa	Signature of technician:						
Next	Next regular inspection:						

Info To print out this document, please use the corresponding PDF file on the CD-ROM.

Help **Contents** 80 of 88

### Spare part order

Page \_\_\_ of \_\_\_

Address:				Volker AG Service Wullener Feld 79 58454 Witten/Gern	nany	LKEF	
Contact perso	on:			Tel.: +49 2302 96096 Fax: +49 2302 96096			
Street addres	s::			E-Mail: service@voell			
Zip code/city/	country:			Shipping address, if	differently from the billing	address	
Telephone nu	mber:		Stamp	Address:			_
Customer nur	mber:						_
Customer ord	ler number:			Attention of:			_
Purchase ord	er date:			Street address:			_
Signature: (Please, to fill in the f	form use block letters)			Zip code/city/country:			_
Please, fill out all in MODEL (Bed-type)	nformation carefully and complete in the form,  SERIAL NOJ YEAR OF MANUFACTURE  (Identifaction label inside of the bed-head)	otherwise we have problems in delivery a	and processing this order.	ITEM P	IUMBER	QUANTITY	
							Info To print out this
							document, please use
							the cor- responding PDF file
							on the CD-ROM.
							_

Contents 81 of 88 Help

### Requirement for using the electronic instructions for use

In order to open the electronic version of these operating instructions (PDF file) found on the CD-ROM, you must have Adobe Reader 7.0.5 or later (or the corresponding Adobe Acrobat version) installed on your PC.

Adobe Reader is available for nearly all operating systems. The newest version can be obtained free of charge by download from <a href="https://www.adobe.de/products/acrobat/readstep2.html">www.adobe.de/products/acrobat/readstep2.html</a> \$\square\$ .

CD-ROM

#### **Trademarks**

Adobe, Acrobat and Reader are protected or registered trademarks of Adobe Systems Incorporated.

Windows is a registered trademark of the Microsoft Corporation.

Notices:

Notices:

Notices:

## **UULKER**

### Völker AG

Wullener Feld 79 · 58454 Witten/Germany Telephone +49 2302 960 96-0 Fax +49 2302 960 96-16 www.voelker.de · info@voelker.de