



# Bedside cabinets and server for hospitals and healthcare facilities

# Völker bedside cabinets

## The ideal complement to any Völker bed.

In addition to the bed, the bedside cabinet is a key aid for patients and residents during their stay. On the one hand, it ensures that all important items, such as books, glasses and electronic devices, are within reach at all times. On the other hand, it enables patients to eat meals in a comfortable and convenient position. Here you will find an overview of the Völker bedside cabinet models.



Model K 940



1× bottle holder



3× bottle holder



Bedpan holder



Attachment point for telephone cradle



Towel holder



3-sided continuous railing



Model	K 931	K 940	K 941	K 941/2	K 941 MK	K941 4.0
Dimensions in cm (W × D × H)	48×48×73	60×48×56	60×48×73	60×48×85	68×56×79	68×48×73
Overbed tray in cm (W × D)	-	55×35	55×35	55×35	55×35	55×35
Height adjustment range overbed tray in cm	-	64-98	64-98	64-98	64-98	64-98
Weight in kg	21	29	29	32	48	35
Material	RAL 9006 / HPL *	RAL 9006 / HPL *	RAL 9006 / HPL *	RAL 9006 / HPL *	RAL 9006 / HPL *	RAL 9006 / HPL *

## Configuration

Bedside table	-	✓	✓	✓	✓	○
Can be used on both sides	✓	✓	✓	✓	-	✓
1× bottle holder, door (on both sides)	○	○	○	○	-	○
3× bottle holder, drawer (on both sides)	○	○	○	○	-	○
Castors Ø 6.5 cm (2 castors with brake)	✓	✓	✓	✓	✓	✓
Door	✓	✓	✓	✓	-	✓
Drawer instead of door	○	○	○	○	-	○
Open storage shelf	✓	-	✓	✓	-	✓
Bedpan holder	○	○	○	○	○	-
Two attachment points for telephone cradle	○	○	○	○	○	○
Tabletop extender	-	✓	-	-	-	-
1 towel holder	○	○	○	○	○	-
Top drawer, lockable	○	○	○	○	○	○
3-sided continuous railing, height: 68 mm 3D	○	-	○	○	○	-



Model	P 920	P 930	P 935	P 941	P 941 O.T.	P 941 4.0
Dimensions in cm (W × D × H)	42×44×61	48×48×61	48×48×61	60×48×74	48×48×74	60×48×73
Overbed tray in cm (W × D)	-	-	-	55×35	-	55×35
Height adjustment range overbed tray in cm	-	-	-	64-98	-	64-98
Weight in kg	21	23	24	32	24	35
Material	Decorative chipboard*	Decorative chipboard*	Decorative chipboard*	Decorative chipboard*	Decorative chipboard*	Decorative chipboard*
Design	-	Round profile	Square profile	Round profile	Round profile	-

## Configuration

Bedside table	-	-	-	✓	-	○
Can be used on both sides	-	-	-	✓	✓	✓
1× bottle holder, door (on both sides)	-	-	-	○	○	○
1× bottle holder on door	○	○	○	-	-	-
3× bottle holder in drawer	○	○	○	-	-	○
Castors Ø 6.5 cm (2 castors with brake)	✓	✓	✓	✓	✓	✓
Drawer instead of door	✓	✓	✓	-	-	-
Open storage shelf	-	-	-	✓	✓	✓
Top drawer, lockable	○	○	○	○	○	○

# Server

## Model

RT-985

Dimensions in cm (W × D × H), without railing	98 × 45 × 74 up to 117
Weight in kg	20
Material	Aluminium (RAL 9006) and high-pressure laminate (HPL)
Colour	See colour collection

## Configuration

Castors Ø 5.0 cm *	×
Castors Ø 6.5 cm *	✓
Castors Ø 7.5 cm *	○
Full-length tray	×

✓ Standard   ○ optional   × optional at no additional cost   \* 2 castors with brake



24/7

**ALWAYS  
AVAILABLE**



**FOCUS  
ON SERVICE**



**HIGH-QUALITY  
MATERIALS**



**CERTIFIED  
SUSTAINABILITY**

**Head office**

**Völker GmbH**

Wullener Feld 79  
58454 Witten  
Germany  
+49 (0) 2302 960 960

**Völker Schweiz GmbH**

Moosstrasse 2  
6003 Luzern  
Switzerland  
+41 (0) 44 563 1939

**Völker Medizinprodukte GmbH**

Kaufmannstraße 14  
6020 Innsbruck  
Austria  
+43 (0) 1 205 8365

**Voelker France SAS**

27 Rue Maurice Flandin  
69444 Lyon Cedex 03  
France  
+33 (0) 18 288 1530

**Voelker Nederland B.V.**

Laan 21  
1741 EA Schagen  
Netherlands  
+31 (0) 20 808 1603

**General enquiries**

info@voelker.de  
+49 (0) 2302 960 960

**Sales / Service / Replacement parts**

kundenservice@voelker.de  
+49 (0) 2302 960 9630