



Bedside cabinets and server for hospitals and healthcare facilities

Völker bedside cabinets

The ideal complement to any bed.

In addition to the bed, the bedside cabinet is a key aid for patients and residents during their stay. On the one hand, it ensures that all important items, such as books, glasses and electronic devices, are within reach at all times. On the other hand, it enables patients to eat meals in a comfortable and convenient position. Here you will find an overview of the Völker bedside cabinet models.

Practical solutions for nursing, mealtimes or relaxation: our bedside cabinets with their optional, height-adjustable bedside tables or the Server RT-985 are extremely versatile and convenient.



1× bottle holder



3× bottle holder



Bedpan holder



Attachment point for telephone cradle



Towel holder



3-sided continuous railing



Hospital Models	K 931	K 940	K 941	K 941/2	K 941 MK
Dimensions in cm (W × D × H)	48×48×73	60×48×73	60×48×73	60×48×86	68×56×79
Overbed tray in cm (W × D)	-	55×35	55×35	55×35	55×35
Height adjustment range overbed tray in cm	-	63-95	63-95	63-95	63-95
Weight in kg	25	35	35	38	48
Material	RAL 9006 / HPL *	RAL 9006 / HPL *	RAL 9006 / HPL *	RAL 9006 / HPL *	RAL 9006 / HPL *

Configuration

Bedside table	-	✓	✓	✓	✓
Can be used on both sides	✓	✓	✓	✓	-
1× bottle holder, door (on both sides)	○	○	○	○	-
3× bottle holder, drawer (on both sides)	○	○	○	○	-
Castors Ø 6.5 cm (2 castors with brake)	✓	✓	✓	✓	✓
Door	✓	✓	✓	✓	-
Drawer instead of door	○	○	○	○	-
Open storage shelf	✓	-	✓	✓	-
Bedpan holder	○	○	○	○	○
Two attachment points for telephone cradle	○	○	○	○	○
Tabletop extender	-	✓	-	-	-
1 towel holder	○	○	○	○	○
Top drawer, lockable	○	○	○	○	○
3-sided continuous railing, height: 68 mm 3D	○	-	○	○	○



Healthcare Models

	P 920	P 930	P 941	P 941 O.T.
Dimensions in cm (W × D × H)	42×44×61	48×48×61	60×48×74	48×48×74
Overbed tray in cm (W × D)	-	-	55×35	-
Height adjustment range overbed tray in cm	-	-	63-95	-
Weight in kg	24	24	33	28
Material	Decorative chipboard*	Decorative chipboard*	Decorative chipboard*	Decorative chipboard*
Design	-	Round profile	Round profile	Round profile

Configuration

Bedside table	-	-	✓	-
Can be used on both sides	-	-	✓	✓
1× bottle holder, door (on both sides)	-	-	○	○
1× bottle holder on door	○	○	-	-
3× bottle holder in drawer	○	○	-	-
Castors Ø 6.5 cm (2 castors with brake)	✓	✓	✓	✓
Drawer	✓	✓	-	-
Door	○	○	✓	✓
Open storage shelf	-	-	✓	✓
Top drawer, lockable	○	○	○	○

Server

Model

RT-985

Dimensions in cm (W × D × H), without railing	98 × 45 × 74 up to 117
Weight in kg	20
Material	Aluminium (RAL 9006) and high-pressure laminate (HPL)
Colour	See colour collection

Configuration

Castors Ø 5.0 cm *	×
Castors Ø 6.5 cm *	✓
Castors Ø 7.5 cm *	○
Full-length tray	○

✓ Standard ○ optional × optional at no additional cost

* 2 castors with brake



Our Server RT-985, a height-adjustable table on wheels is a great all-rounder – position it over the bed, beside the bed, or anywhere else in the room. For hospital and healthcare.

24/7

**ALWAYS
AVAILABLE**



**FOCUS
ON SERVICE**



**HIGH-QUALITY
MATERIALS**



**CERTIFIED
SUSTAINABILITY**

Head office

Völker GmbH

Wullener Feld 79
58454 Witten
Germany
+49 (0) 2302 960 960

Völker Schweiz GmbH

Moosstrasse 2
6003 Luzern
Switzerland
+41 (0) 44 563 1939

Völker Medizinprodukte GmbH

Kaufmannstraße 14
6020 Innsbruck
Austria
+43 (0) 1 205 8365

Voelker France SAS

27 Rue Maurice Flandin
69444 Lyon Cedex 03
France
+33 (0) 18 288 1530

Voelker Nederland B.V.

Laan 21
1741 EA Schagen
Netherlands
+31 (0) 20 808 1603

General enquiries

info@voelker.de
+49 (0) 2302 960 960

Sales / Service / Replacement parts

kundenservice@voelker.de
+49 (0) 2302 960 9630