



# Healthcare bed accessories

2080 / 3080 / 5384

# Völker accessories

## Healthcare beds 2080, 3080 and 5384.

All Völker accessories can be attached by simply inserting or hooking them into the integrated mounting sleeves, brackets and rails.

The following pages provide a detailed overview of the accessories available for each healthcare bed.

## Three convincing elements

- ✓ Völker offers a wide range of practical accessories that are perfectly tailored to Völker beds.
- ✓ All accessories are designed to be easily inserted or hooked into integrated universal rails.
- ✓ All healthcare beds feature a universal rail, and as such the side rails never prevent the usage of accessories.

## Völker healthcare beds - outstanding benefits



**SAFE  
WORKINGLOAD**



**WIDE  
ADJUSTMENT RANGE**




**DYNAMIC  
EXTENSION**





**LOW  
BED HEIGHT**



**LOW POSITION  
FOR FEET**



People requiring care have a particular right to respect and regard for their dignity. Caregivers also have a right to equipment that makes their daily work as easy as possible. We think, develop and manufacture with this in mind.

Configuration		Item no.	2080	3080	5384
	<b>Trapeze bar with handle</b> Material: powder coated	PZP-3060	✓	✓	✓
	<b>Trapeze bar with handle</b> Adjustable with one hand Material: powder coated	PZP-3060/2	✓	✓	✓
	<b>Assist rail spacer</b> 1 piece Material: Stainless steel and high-pressure laminate (HPL)	PZK-969/2	-	✓	✓
	<b>Assist rail spacer cover</b> 1 piece Colour: Brown	PZP-3059	-	✓	✓
	<b>Assist rail spacer holder</b>	PZP-3070	-	✓	✓
	<b>Fixation belt holder</b> 1 pair				
	Relocatable without tools	PZK-972/3	✓	✓	-
	Fixed in place by screws	PZK-972/4	✓	✓	-
	<b>Crutch holder</b> 1 pair	PZP-965	✓	✓	✓
	<b>Hand control holder</b>	PZK-1044/3	-	✓	-
	<b>Care Clip cable holder</b> For cable routing on the trapeze bar e.g. for a nurse call button	KT3860	✓	✓	✓
	<b>Infusion holder</b> For trapeze bar PZP-3060	PZP-3055/2	✓	✓	✓
	<b>Infusion stand</b> Adjustable with one hand 2 kg load capacity per hook Material: stainless steel *PZP-5084 adapter required	PZP-3057	✓	✓	✓
	<b>ISO standard rail</b> Material: stainless steel				
	Length: 40 cm	PZK-1083/6	✓	✓	-
	Length: 60 cm	PZK-1083/8	✓	✓	-

✓ Standard - not available \* not with fixed side rails

Configuration		Item no.	2080	3080	5384
	<b>Side rail protective cover</b> For bed length: 200 cm Height of the side rail: 34 cm Colour: Brown	PZP-2056-2.0	✓	-	-
	<b>Side rail protective covers</b> For bed length: 200 cm Height of the side rail: 34 cm Colour: Brown	PZP-3056_C_2.0	-	✓	✓
	<b>Universal hook</b> 1 pair	PZK-1043	✓	✓	-
	<b>Universal hook</b> 1 pair For ISO standard rails	PZK-943	✓	✓	✓
	<b>Urine bottle basket</b> For hanging on the accessories rail Material: Chrome	PZK-1040	✓	✓	-
	<b>Urine bottle basket</b> Inserts into the mounting sleeve Material: Chrome	PZK-1040/2	-	✓	✓
	<b>Wall protector wheels</b> Vertical	PZP-3068	✓	✓	✓
	<b>Maya lamp</b> Connected behind the headboard using an adapter	PMAYA02	✓	✓	✓
	<b>Tablet holder</b> Can be plugged into the socket of the hand control	PZK-1092	-	✓	✓
	<b>Mobile phone holder</b> Can be plugged into the socket of the hand control	PZK-1093	-	✓	✓

With some accessories, availability also depends on the year of manufacture and configuration of your bed. Need to know exactly which accessory versions are suitable for your Völker bed? Then please check the serial number of your bed and contact your sales partner.






✓ Standard - not available \* not with fixed side rails

# All Völker healthcare beds at a glance

Technical data.








2080

Safe working load		230 kg
Unladen weight		approx. 96 kg
Lying surface width		90 cm
Lying surface length		200 cm
External dimensions		98 cm








3080

Safe working load		230 kg
Unladen weight		approx. 96 kg
Lying surface width		90 / 100 cm
Lying surface length		200 / 210 cm
External dimensions		98 cm



5384

Safe working load		230 kg
Unladen weight		127 – 158 kg
Lying surface width		90 cm
Lying surface length		200 cm
External dimensions		98 cm

Learn more about  
the technical data  
[voelker.de](http://voelker.de)

24/7

**ALWAYS  
AVAILABLE**



**FOCUS  
ON SERVICE**



**HIGH-QUALITY  
MATERIALS**



**CERTIFIED  
SUSTAINABILITY**

**Head office**

**Völker GmbH**

Wullener Feld 79  
58454 Witten  
Germany  
+49 (0) 2302 960 960

**Völker Schweiz GmbH**

Moosstrasse 2  
6003 Luzern  
Switzerland  
+41 (0) 44 563 1939

**Völker Medizinprodukte GmbH**

Kaufmannstraße 14  
6020 Innsbruck  
Austria  
+43 (0) 1 205 8365

**Voelker France SAS**

27 Rue Maurice Flandin  
69444 Lyon Cedex 03  
France  
+33 (0) 18 288 1530

**Voelker Nederland B.V.**

Laan 21  
1741 EA Schagen  
Netherlands  
+31 (0) 20 808 1603

**General enquiries**

info@voelker.de  
+49 (0) 2302 960 960

**Sales / Service / Replacement parts**

kundenservice@voelker.de  
+49 (0) 2302 960 9630