



Hospital bed

S 962-2 Vis-a-Vis

Völker S 962-2 Vis-a-Vis

Mobilisation made easy.

Not every patient who cannot or does not want to leave their bed is completely immobile as a result. Often, physical ailments or medical interventions prevent patients from getting out of bed via a lateral transfer. Yet mobilisation is also possible if the bed can be adjusted so that it forms a sitting position that is open to the front and the patient's feet are immediately in stable contact with the floor. This means that – depending on the patient's resources – the process of regaining mobility can occur in stages. From lying to sitting to standing and walking. Through a successful combination of form and function, Völker hospital beds provide greater patient comfort and enable care staff to work safely and efficiently. Völker's unique healthcare design ensures that beds fulfil their medical purpose and at the same time create a homelike ambience that supports patients' well-being and recovery.

Four convincing elements

- ✓ Forward open seating position for mobility and care at eye level.
- ✓ Split, height-adjustable side rails – move along with the lying surface as it is adjusted.
- ✓ With the innovative Völker MiS® Micro-Stimulation System. For greater patient comfort and better sleep.
- ✓ Relief for the care staff through independent mobilisation for the patient.



**SAFE
WORKINGLOAD**



**WIDE
ADJUSTMENT RANGE**



**ERGONOMIC
WORKING HEIGHT**



**LOW POSITION
FOR FEET**








The sitting position can be reached effortlessly, so that patients can talk to care staff and visitors face to face, at eye level.

Safe mobilisation for patients

Relief for the care staff thanks to the patented Völker side rail.

Independent mobilisation for patients helps with recovery following a period of illness. This mobilisation is supported by the central gap in the Völker side rail, alleviating the restricted freedom of movement imposed by side guards. The stable, multi-telescopic design of the Völker side rail enables patients to stand up independently and regain confidence in their own bodies. As standard, the Völker side rails are fitted on all sides of the lying surface and offer optimum protection in every position. A fixed split version of the side rail is also available as an option.

Technical data

Safe working load		250 kg
Unladen weight		approx. 169 kg
Lying surface width		90 cm
Lying surface length		210 cm
External dimensions		99 cm



Learn more about
the technical data
voelker.de



Greater comfort and relaxation

The micro-stimulation lying surface.

All of Völker's lying surfaces feature the patented MiS® Micro-Stimulation System, which is based on neuroscientific principles. It reacts to even the most minute of the patient's own movements, relaying them back to the body as tactile stimuli from the lying surface through the mattress. Experience in recent years has shown that such feedback can lead to improved perception of the body, as well as improvements in sleeping and lying comfort. The combination of lying surface and mattress may also reduce the risk of bedsores.¹



Völker offers a wide range of practical accessories that are perfectly tailored to Völker beds.



Different head- and footboards available.

Finess mattress

Two-section construction for Vis-a-Vis beds.

- ✓ Highly elastic HR cold foam mattress with 3-layer design
- ✓ Firm elasticated, stabilising middle layer
- ✓ Foam materials: certified according to Oeko-Tex Standard 100 product class 1
- ✓ Special wave cut for effective temperature regulation
- ✓ Reduction in moisture build up and promotion of air circulation through additional air channels
- ✓ The mattress can be turned in four directions due to the symmetrical design
- ✓ The mattress with cover can be disinfected in an autoclave at temperatures up to 105 °C

Configuration

	Finess
Colour: grey	✓
Textile carrier: 100% polyester, Wetness barrier layer: 100% polyurethane	✓
Bi-elastic	✓
Bacteria, mite and virus-proof	✓
Suitable for allergy sufferers	✓
Anti-bacterial coating: Preventol MP100 with Lanxess	✓
Water and urine impermeable	✓
Flame retardant in accordance with UNI 9175: Class 1M + BS 7175 crib5+7	✓
Resistant to blood, urine, oils, fats	✓
Resistant to disinfectants according to VAH & RKI	✓
Washable up to 95 °C	✓
Autoclavable up to 105 °C	✓
Surface weight: 190 g/m ² ± 10	✓
Water column according to UNI EN 20811 > 200 cm	✓
Water vapour permeable according to DIN 53122 > 500 g/m ² /24 hour	✓
Continuous zip over one long side and one narrow side	✓
Removable cover	✓



24/7

**ALWAYS
AVAILABLE**



**FOCUS
ON SERVICE**



**HIGH-QUALITY
MATERIALS**



**CERTIFIED
SUSTAINABILITY**

Head office

Völker GmbH

Wullener Feld 79
58454 Witten
Germany
+49 (0) 2302 960 960

Völker Schweiz GmbH

Moosstrasse 2
6003 Luzern
Switzerland
+41 (0) 44 563 1939

Völker Medizinprodukte GmbH

Kaufmannstraße 14
6020 Innsbruck
Austria
+43 (0) 1 205 8365

Voelker France SAS

27 Rue Maurice Flandin
69444 Lyon Cedex 03
France
+33 (0) 18 288 1530

Voelker Nederland B.V.

Laan 21
1741 EA Schagen
Netherlands
+31 (0) 20 808 1603

General enquiries

info@voelker.de
+49 (0) 2302 960 960

Sales / Service / Replacement parts

kundenservice@voelker.de
+49 (0) 2302 960 9630