



# Hospital bed

## S 966

# Völker S 966

## One bed. More options.

Meeting the needs of patients and simplifying the daily work of care staff – the Völker S 966 hospital bed fulfils both of these objectives in equal measure.

The Völker S 966 provides maximum safety for patients. In addition, the unique Völker healthcare design promotes a homelike atmosphere, thereby supporting a patient's well-being and recovery.

With the Völker S 966, care staff benefit from numerous improvements that simplify handling and save time. This brings noticeable relief in everyday hospital life.

The Völker S 966 is designed in such a way that it can easily cope with all typical hospital situations. Your hospital will benefit from the flexible utilisation together with the consistent and attractive appearance.

## Five convincing elements

- ✓ Meets the highest standards of quality, functionality and safety.
- ✓ Available in numerous versions – modular elements are easy to combine.
- ✓ With numerous innovative Völker product functions as standard.
- ✓ Successful synthesis of form and function.
- ✓ Optionally available with a central fifth castor.



**SAFE  
WORKINGLOAD**



**WIDE  
ADJUSTMENT RANGE**




**DYNAMIC  
EXTENSION**



**LOW  
BED HEIGHT**



**TRENDELENBURG**



For increased well-being and a faster recovery: The Völker S 966 is based on Völker's modern, comfort-oriented healthcare design. High-quality materials create a pleasant look and feel.








# Safe mobilisation for patients

## Relief for the care staff thanks to the patented Völker side rail.

Independent mobilisation for patients helps with recovery following a period of illness. This mobilisation is supported by the central gap in the Völker side rail, alleviating the restricted freedom of movement imposed by side guards. The stable, multi-telescopic design of the Völker side rail enables patients to stand up independently and regain confidence in their own bodies. As standard, the Völker side rails are fitted on all sides of the lying surface and offer optimum protection in every position. A fixed split version of the side rail is also available as an option.

### Technical data

Safe working load		270 kg
Unladen weight		approx. 160 kg
Lying surface width		90 / 100 cm
Lying surface length		200 / 210 cm
External dimensions		99 – 104 cm*

\*for fixed split side rails




#### Central brake

Can also be operated at the lowest bed position. This means that the bed can be moved at any height setting.



Learn more about  
the technical data  
[voelker.de](https://voelker.de)



## Greater comfort and relaxation

### The micro-stimulation lying surface.

All of Völker's lying surfaces feature the patented MiS® Micro-Stimulation System, which is based on neuroscientific principles. It reacts to even the most minute of the patient's own movements, relaying them back to the body as tactile stimuli from the lying surface through the mattress. Experience in recent years has shown that such feedback can lead to improved perception of the body, as well as improvements in sleeping and lying comfort. The combination of lying surface and mattress may also reduce the risk of bedsores.<sup>1</sup>



#### Backlit hand control

Can be easily hung off the side rail and is easy to move.



#### Automatic features

Defined positions can be set quickly and easily with a simple double click. Enables users to better focus on care activities.





# The better choice

For patients and care staff.

- ✓ Low bed height of < 30 cm.  
Significantly reduces the severity of injuries in the event of a fall.
- ✓ Head- and footboards can be removed quickly and without tools, in an emergency, for example, or for cleaning.
- ✓ Large selection of colours and decors to suit your requirements.



Völker offers a wide range of practical accessories that are perfectly tailored to Völker beds.



People requiring care have a particular right to respect and regard for their dignity. Caregivers also have a right to equipment that makes their daily work as easy as possible. We think, develop and manufacture with this in mind.

24/7

**ALWAYS  
AVAILABLE**



**FOCUS  
ON SERVICE**



**HIGH-QUALITY  
MATERIALS**



**CERTIFIED  
SUSTAINABILITY**

**Head office**

**Völker GmbH**

Wullener Feld 79  
58454 Witten  
Germany  
+49 (0) 2302 960 960

**Völker Schweiz GmbH**

Moosstrasse 2  
6003 Luzern  
Switzerland  
+41 (0) 44 563 1939

**Völker Medizinprodukte GmbH**

Kaufmannstraße 14  
6020 Innsbruck  
Austria  
+43 (0) 1 205 8365

**Voelker France SAS**

27 Rue Maurice Flandin  
69444 Lyon Cedex 03  
France  
+33 (0) 18 288 1530

**Voelker Nederland B.V.**

Laan 21  
1741 EA Schagen  
Netherlands  
+31 (0) 20 808 1603

**General enquiries**

info@voelker.de  
+49 (0) 2302 960 960

**Sales / Service / Replacement parts**

kundenservice@voelker.de  
+49 (0) 2302 960 9630