



Hospital bed accessories

S 962-2 / S 962-2 W

S 962-2 Vis-a-vis / S 966

Völker accessories

for hospital beds S 962-2, S 962-2 W,
S 962-2 Vis-a-Vis and S 966.

All Völker accessories can be attached by simply inserting or hooking them into the integrated mounting sleeves, brackets and rails.

The following pages provide a detailed overview of which accessories are available for which hospital beds.

Three convincing elements

- ✓ Völker offers a wide range of practical accessories that are perfectly tailored to Völker beds.
- ✓ All accessories are designed to be easily inserted or hooked into integrated universal rails.
- ✓ All hospital beds feature a universal rail, and as such the side rails never prevent the usage of accessories.

Völker hospital beds – Outstanding benefits



**SAFE
WORKINGLOAD**



**WIDE
ADJUSTMENT RANGE**



**DYNAMIC
EXTENSION**













**LOW
BED HEIGHT**









People requiring care have a particular right to respect and regard for their dignity. Caregivers also have a right to equipment that makes their daily work as easy as possible. We think, develop and manufacture with this in mind.

Configuration		Item no.	S 962-2 S 962-2 W	S 962-2 VaV	S 966
	Fold-down head- or footboard Only available with a bed width of 90 cm	PZK-1035/4	✓*	-	✓*
	Trapeze bar with handle Material: powder coated	PZP-3060	✓	-	✓
	Trapeze bar with handle Material: stainless steel	PZK-936/1	✓	✓	✓
	Trapeze bar with handle Adjustable with one hand				
	Material: powder coated	PZP-3060/2	✓	-	✓
	Material: stainless steel	PZK-936/3	✓	✓	✓
	Document holder	PZK-1070	✓	✓	✓
	Assist rail spacer 1 piece Material: Stainless steel and high-pressure laminate (HPL)	PZK-969/2	✓*	✓	✓*
	Assist rail spacer cover Colour: grey	PZK-982	✓*	✓	✓*
	Assist rail spacer holder	PZK-979	✓*	✓	✓*
	Universal belt holder 1 pair				
	Relocatable without tools	PZK-972/3	✓	✓	✓
	Fixed in place by screws	PZK-972/4	✓	✓	✓

Configuration		Item no.	S 962-2 S 962-2 W	S 962-2 VaV	S 966
	Crutch holder 1 pair	PZK-965	✓	✓	✓
	Crutch holder 1 holder fits 2 crutches	PZK-965/2	✓	✓	✓
	Hand control holder	PZK-1044/3	✓*	✓	✓*
	Care Clip cable holder for cable routing on the trapeze bar				
	For trapeze bar PZP-936 (38 mm)	KT3860	✓	✓	✓
	For trapeze bar PZK-3060 (31 mm)	KT3895	✓	✓	✓
	Infusion holder For trapeze bar PZP-3060	PZP-3055/2	✓	-	✓
	Infusion holder For trapeze bar PZK-936	PZK-937/3	✓	✓	✓
	Infusion stand Adjustable with one hand Load capacity per hook: 2 kg Material: stainless steel	PZK-939_R	✓	✓	✓
	Infusion stand Adjustable with one hand Load capacity per hook: 2 kg Material: stainless steel	PZK-939_S	✓	✓	✓
	Infusion stand (straight) Adjustable with one hand Load capacity per hook: 2 kg Material: stainless steel	PZK-939_U	✓	✓	✓
	ISO standard rail Material: stainless steel				
	Length: 40 cm	PZK-1083/6	✓	✓	✓
	Length: 60 cm	PZK-1083/8	✓	✓	✓

✓ available - not available * not with fixed side rails

Configuration		Item no.	S 962-2 S 962-2 W	S 962-2 VaV	S 966
	ISO standard rail Fixed to the backrest at the headboard Material: stainless steel	PZK-1091	✓	✓	-
	Name plate	PZK-950_R	✓	✓	✓
	Equipotential cable Length: 300 cm	PZK-948	✓	✓	✓
	Redon drain basket Material: Chrome	PZK-1041	✓*	✓	-
	Writing surface Material: High-pressure laminate (HPL) Colour: grey	PZK-962/2	✓	✓	✓
	Towel holder	PZK-978/2	✓	✓	✓
	Side rail guards For bed length: 200 cm				
	Height of the side rail: 34 cm	PZK-981_2.0	✓*	✓	✓*
	Height of the side rail: 37 cm	PZK-981EN20	✓*	✓	✓*
	For fixed side rails	PZK-984	✓	-	✓
	Universal hook 1 pair	PZK-1043	✓	✓	✓
	Universal hook 1 pair For ISO standard rails: PZK-1083/PZK-1064	PZK-943	✓	✓	✓

Configuration		Item no.	S 962-2 S 962-2 W	S 962-2 VaV	S 966
	Urine bottle basket For hanging on the accessories rail Material: Chrome	PZK-1040	✓*	✓	-
	Urine bottle basket Inserts into the mounting sleeve Material: Chrome	PZK-1040/2	✓*	✓	✓*
	Urine bottle basket For ISO standard rails: PZK-1083 Material: Chrome	PZK-940/3	✓	✓	-
	Maya lamp Connected behind the headboard using an adapter	PMAYA01	✓	✓	✓
	Tablet holder Can be plugged into the socket of the hand control	PZK-1092	✓	✓	✓
	Mobile phone holder Can be plugged into the socket of the hand control	PZK-1093	✓	✓	✓
Arm rests 1 set		PZK-ARVAV	-	✓	-
Foot rests 1 set		PZK-FRVAV	-	✓	-
Bed sheets Two-part sheets for all bed sizes		PZK-BSVAV	-	✓	-

With some accessories, availability also depends on the year of manufacture and configuration of your bed. Need to know exactly which accessory versions are suitable for your Völker bed? Then please check the serial number of your bed and contact your sales partner.

All Völker hospital beds at a glance

Technical data.



S 962-2 / S 962-2 W

Safe working load		250 kg
Unladen weight		approx. 143 / 164 kg
Lying surface width		90 / 100 cm
Lying surface length		200 / 210 cm
External dimensions		99 - 104 cm*



S 962-2 Vis-a-Vis

Safe working load		250 kg
Unladen weight		approx. 169 kg
Lying surface width		90 cm
Lying surface length		210 cm
External dimensions		99 cm



S 966

Safe working load		270 kg
Unladen weight		approx. 160 kg
Lying surface width		90 / 100 cm
Lying surface length		200 / 210 cm
External dimensions		99 - 104 cm*

*for fixed split side rails

Learn more about
the technical data
voelker.de

24/7

**ALWAYS
AVAILABLE**



**FOCUS
ON SERVICE**



**HIGH-QUALITY
MATERIALS**



**CERTIFIED
SUSTAINABILITY**

Head office

Völker GmbH

Wullener Feld 79
58454 Witten
Germany
+49 (0) 2302 960 960

Völker Schweiz GmbH

Moosstrasse 2
6003 Luzern
Switzerland
+41 (0) 44 563 1939

Völker Medizinprodukte GmbH

Kaufmannstraße 14
6020 Innsbruck
Austria
+43 (0) 1 205 8365

Voelker France SAS

27 Rue Maurice Flandin
69444 Lyon Cedex 03
France
+33 (0) 18 288 1530

Voelker Nederland B.V.

Laan 21
1741 EA Schagen
Netherlands
+31 (0) 20 808 1603

General enquiries

info@voelker.de
+49 (0) 2302 960 960

Sales / Service / Replacement parts

kundenservice@voelker.de
+49 (0) 2302 960 9630