



Hospital beds

Protect

Völker. Caring for People

VÖLKER

No compromise on safety

Tailor-made solutions for special applications.

The safety of both patients and caregivers is Völker's top priority. With this in mind, we have optimised the existing S 962-2 and S 966 hospital beds specifically for areas where patients need a higher level of protection or security.

The features of Völker beds enable caregivers to carry out their tasks quickly and easily, allowing them to focus fully on the patient rather than on how to operate the bed. The beds differ only in their standard configuration. For secure areas, the beds are available with the same safety features – either as a complete package or configured with individual features.

Völker hospital beds – Outstanding benefits



**SAFE
WORKINGLOAD**



**WIDE
ADJUSTMENT RANGE**



**DYNAMIC
EXTENSION**



**LOW
BED HEIGHT**

Our Protect collection

Key facts at a glance.



S 962-2 Protect

Safe working load		250 kg
Unladen weight		approx 143 kg
Lying surface width		90 cm
Lying surface length		200 cm
External dimensions		99 cm
Height		approx 40 – 80 cm*



S 966 Protect

Safe working load		270 kg
Unladen weight		approx 160 kg
Lying surface width		90 cm
Lying surface length		200 cm
External dimensions		99 cm
Height		approx 28,5 – 88,5 cm*

* depending on configuration



Safe and reliable

Protect configuration features.

In order to ensure greater safety for patients and caregivers, the S 962-2 and S 966 hospital beds have been optimised and now include additional features.

To minimise the sources of danger and the potential misuse of the bed by the patient, the head- and footboards as well as the HPL lying surface are fixed in place by screws. The hand controls can also be removed and the nurse keypads can be locked, such that both are out of patients' reach.



Power cable in Anaconda cable housing

Suicidal patients need a greater level of protection. For this reason, the power cable is enclosed in Anaconda cable housing, which reduces the possibility of misuse.

Simple operation means greater focus on patient care

The enhanced version of the Völker classic S 962-2.

Easy operation of the bed is particularly important in areas where patients require a higher level of protection. The low weight of the S 962-2 Protect means that it can be easily moved, saving the nursing staff both time and stress. In addition, the four lifting columns make the S 962-2 Protect stable and simple to clean, and allow it to be lowered easily.

More solutions for more safety

The S 966 low-height bed reduces the risk of falling.

The low-height bed S 966 has been designed to reduce injuries caused by patients falling out of bed, specifically in areas where patients require a higher level of protection. Depending on the configuration, the S 966 Protect can be lowered to a height of approx. 28.5 cm and raised to an ergonomic working height of 88.5 cm. A fifth castor can be fitted to the bed if desired.



Universal belt holder

Restraining a patient is usually the last resort to protect both the patient and caregivers. Twelve die-cast zinc universal belt holders are screwed onto the side of the bed so that the patient can be restrained quickly and securely.

Configuration

	S962-2 Protect	S966 Protect
Design		
Design S	✓	✓
Colours		
in accordance with the Völker colour collection	✓	✓
Lying surface width		
90 cm	✓	✓
Lying surface length		
200 cm	✓	✓
Lying surfaces		
HPL (screwed on)	✓	✓
Height of side rails		
Moving split side rail: 34- 35.5 cm	✓	✓
Moving split side rail: 37- 38.5 cm	○	○
Moving split side rail: 43.5-45 cm	○	○
Mounting fixtures		
Single mount - external	○	○
Controls		
Hand control with hook, longitudinal, attachable both sides	✓	✓
Hand control with hook, transverse, attachable both sides	○	○
Nurse keypad, can be locked away in the linen holder	○*	○*
Castors		
Single castors S Ø125 mm	-	✓
Single castors S Ø150 mm	✓	○
Integral single castors Ø 150 mm	○	○
Linea twin castors Ø125 mm	-	○
Twin castors Ø150 mm	○	○
Linea twin castors Ø150 mm	-	○
Track caster		
at the foot of the bed	✓	✓
at the head of the bed	○	○
5. Castor (only with Linea twin castors Ø150 mm)	-	○
Central brake		
Foot pedal at foot of bed - on both sides	-	✓
Additional configuration options		
Anaconda cable 900 mm	○*	○*
12 die-cast zinc universal belt holder screwed to the lying surface	○*	○*
Head- and footboards fixed in place by screws	○*	○*
Bed extension	○	○
Underbed light	-	○

✓ Standard ○ optional - not available * Standard with complete package

24/7

**ALWAYS
AVAILABLE**



**FOCUS
ON SERVICE**



**HIGH-QUALITY
MATERIALS**



**CERTIFIED
SUSTAINABILITY**

Head office

Völker GmbH

Wullener Feld 79
58454 Witten
Germany
+49 (0) 2302 960 960

Völker Schweiz GmbH

Moosstrasse 2
6003 Luzern
Switzerland
+41 (0) 44 563 1939

Völker Medizinprodukte GmbH

Kaufmannstraße 14
6020 Innsbruck
Austria
+43 (0) 1 205 8365

Voelker France SAS

27 Rue Maurice Flandin
69444 Lyon Cedex 03
France
+33 (0) 18 288 1530

Voelker Nederland B.V.

Laan 21
1741 EA Schagen
Netherlands
+31 (0) 20 808 1603

General enquiries

info@voelker.de
+49 (0) 2302 960 960

Sales / Service / Replacement parts

kundenservice@voelker.de
+49 (0) 2302 960 9630