



# Healthcare bed

## 5384

The 5384 healthcare bed offers maximum protection for residents who are at extreme risk of falling out of bed. Falling out of bed can result in the occupant losing their independence, and can also have other negative effects both on the resident as well as on the healthcare facility. This versatile, low-height bed can be lowered to 24 cm\* and raised to a stable and ergonomic working height of 80 cm to ensure a safe and ergonomic working position for healthcare staff.

- ✓ The moving split side rail ensures sufficient safety in any bed position.
- ✓ The innovative MiS® Micro-Stimulation System provides residents with more comfort and a better sleep.
- ✓ The bed height memorisation feature enables residents and care staff to set the desired height and position with ease.
- ✓ A subtle night light under the bed enhances safety for both residents and caregivers (optional).



## Configuration

5384

<b>Material head- and footboards</b>	
Frame in solid wood, fillings in melamine	✓
<b>Design</b>	
Design S - round profile	✓
Design FS - flat profile (only in combination with 34 - 35.5 cm side rails)	○
<b>Colours</b>	
in accordance with the Völker colour collection	✓
<b>Lying surface width</b>	
90 cm	✓
<b>Lying surface length</b>	
200 cm	✓
<b>Lying surface</b>	
Völker MiS® Micro-Stimulation System <sup>1</sup>	✓
High-pressure laminate	○
<b>Height of the side rail (dependent upon the lying surface)</b>	
Moving split side rail: 34 - 35.5 cm	✓
Moving split side rail: 37 - 38.5 cm	○
Moving split side rail: 43.5 - 45 cm	○
<b>Controls</b>	
Hand control with all function lockout, longitudinal	✓
<b>Castors</b>	
Twin castors, 4 × (Ø 100 mm) - Single castors, 2 × (Ø 80 mm) - head end, lying surface height = 24 - 79 cm	✓
Twin castors, 4 × (Ø 100 mm) - lying surface height = 25 - 81 cm	○
Linea twin castors, 4 × (Ø 125 mm) - lying surface height = 28 - 83 cm	○
Integral S single castors, 4 × (Ø 125 mm) - lying surface height = 31 - 88 cm	○
Integral S single castors, 4 × (Ø 150 mm) - lying surface height = 33 - 88 cm	○
<b>Additional configuration options</b>	
CPR quick lowering	✓
Underbed light	○
Linen holder	○
Bed extension	○
Accessories rail with fixation belt holder on both sides	○
Lateral wall protector wheels, vertically at head end	○

✓ Standard ○ optional \* depending on the selected configuration