



# Hospital bedside cabinet

## K 941 4.0

**Völker bedside cabinets – the ideal complement to any Völker bed. The hospital bedside cabinet offers a number of new optional features and benefits that make the days spent recovering in hospital even more comfortable.**

- ✓ USB ports and sockets (230 V) on both sides
- ✓ Optional ambient lighting
- ✓ Optional underbody lighting with motion sensor to make it easier for patients and care staff to find their way at night
- ✓ Optional bedside table

## Configuration

K 941 4.0

Bedside table	○
Can be used on both sides	✓
1× bottle holder, door, on both sides	○
3× bottle holder, drawer, on both sides	○
Castors 6.5 cm in diameter, 2 castors with brake	✓
Door	✓
Drawer instead of door	○
Open storage shelf	✓
Two possible attachment points for telephone cradle	○
Top drawer, lockable	○
Mains cable	✓
Ambient lighting	○
USB ports on both sides	✓
Power sockets on both sides	✓
Underbody lighting including motion sensor	○

✓ Standard ○ optional



## Technical data

**Width × depth × height** 60 × 48 × 73 cm

**Overbed tray** 55 × 35 cm

**Height adjustment range for the overbed tray** 64 – 98 cm

**Weight** 34 kg

**Material** aluminium (RAL 9006) and high-pressure laminate (HPL)

**Colour** In accordance with the Völker colour collection