

Healthcare bedside cabinet

P 941 4.0

Völker bedside cabinets - the ideal complement to any Völker bed. The new healthcare bedside cabinet offers a number of new optional features and benefits.

- ✓ USB ports on both sides
- ✓ Illuminated drawer (only in combination with ambient lighting)
- ✓ Sockets (230 V) on both sides
- ✓ Underbody lighting with motion sensor to make it easier for residents and care staff to find their way at night
- ✓ Bedside table can be selected as an option

Configuration

P 941 4.0

USB ports on both sides	✓
Power sockets on both sides	✓
Can be used on both sides	✓
Mains cable	✓
Open storage shelf	✓
Door	✓
Castors 6.5cm in diameter, 2 castors with brake	✓
Ambient lighting	○
Underbody lighting including motion sensor	○
Bedside table	○
1× bottle holder, door, on both sides	○
Top drawer, lockable	○

✓ Standard ○ optional



Technical data

Width × depth × height 60 × 48 × 73 cm

Overbed tray 55 × 35 cm

Height adjustment range for the overbed tray 64 - 98 cm

Weight 35 kg

Material decorative chipboard

Colour In accordance with the Völker colour collection