



# Hospital bed

## S 962-2

**The S 962-2 Völker hospital bed meets a wide variety of requirements in hospitals. It enables patients to go through the different stages of mobilisation in a safe way, with support from the care staff – from sitting in the comfort chair position to getting up and leaving the bed.**

- ✓ Two side rail variants: moving split or fixed split side rails (S 962-2 FSG). Both variants are multilevel and can be positioned individually.
- ✓ The innovative MiS® Micro-Stimulation System provides patients with greater comfort and improves the quality of sleep.
- ✓ The comprehensive functionality helps care staff to mobilise patients safely and to reduce the physical strain of their work.
- ✓ Völker's unique healthcare design ensures that hospital beds give the patient's room a homelike atmosphere that promotes well-being and recovery.
- ✓ Also available in version S 962-2 W, which is suitable for washing tunnels.



# VÖLKER

## Configuration

	S 962-2
<b>Design</b>	
Design S	✓
Design S, with fold-down head- or footboard, not available with 100 cm width	○
<b>Colours</b>	
in accordance with the Völker colour collection	✓
<b>Lying surface width</b>	
90 cm	✓
100 cm, not available with fixed headboard	○
<b>Lying surface length</b>	
200 cm	✓
210 cm	○
<b>Lying surface</b>	
Völker MiS® Micro-Stimulation System <sup>1</sup>	✓
High-pressure laminate HPL	○
<b>Height of the side rail (dependent upon the lying surface)</b>	
Moving split side rail: 34 to 35.5 cm	✓
Moving split side rail: 43.5 to 45 cm	○
Moving split side rail: 37 to 38.5 cm	○
Fixed split side rail: 42 to 43.5 cm	○
<b>Mounting fixtures</b>	
Single mount – external	✓
Double mount – internal	○
<b>Controls</b>	
Hand control with hook, longitudinal	✓
Hand control with hook, transverse, not available with fixed split side rails	○
Nurse keypad in linen holder	✓
Nurse keypad in linen holder with single lock with/without automatic functions	○
Nurse keypad in linen holder with enhanced features with/without automatic function	○
<b>Castors</b>	
Integral S single castors (Ø 150 mm)	✓
Integral single castors (Ø 150 mm)	○
Integral twin castors (Ø 150 mm)	○
<b>Track caster</b>	
at the foot of the bed	✓
at the head of the bed	○
<b>Central brake</b>	
Foot pedal at foot end, on both sides	✓
Foot pedal on four sides	○
<b>Additional configuration options</b>	
Bed extension	○
Fixed headboard, not available with 100 cm width and only with external mounting fixture	○

✓ Standard   ○ optional

voelker.de

Völker reserves the right to make changes to the design, specifications and models without prior notice. Products may vary slightly from the colours shown here. Subject to printing errors. All information is non-binding. Not all products/options are available in all countries. For further information on our products or services, please contact your local Völker representative or visit our website.