



Hospital bed

S 966

The S 966 model hospital bed is the result of the continuous further development of the Völker product line. Timeless design combined with innovative technology and Made in Germany.

- ✓ Two side rail variants: moving split or fixed split side rail. Both variants are multi-level and can be positioned individually.
- ✓ The innovative MiS® Micro-Stimulation System provides patients with greater comfort and improves the quality of sleep.
- ✓ Use of a 5th castor (optional)
- ✓ Safe working load 270 kg

Configuration

	S 966
Design S	✓
Colours	
in accordance with the Völker colour collection	✓
Lying surface width	
90 cm	✓
100 cm	○
Lying surface length	
200 cm	✓
210 cm	○
Lying surface	
Völker MiS® Micro-Stimulation System ¹	✓
High-pressure laminate	○
Height of the side rail (dependent upon the lying surface)	
Moving split side rail: 34 – 35.5 cm	✓
Moving split side rail: 37 – 38.5 cm	○
Moving split side rail: 43.5 – 45 cm	○
Fixed split side rail: 42 – 43.5 cm	○
Mounting fixtures	
Double mount – internal	✓
Single mount – external	○
Controls	
Hand control with hook, longitudinal	✓
Hand control with hook, transverse (not with fixed split side rail variant)	○
Nurse keypad in linen holder – with single lock with/without automatic functions	✓
Nurse keypad in linen holder – with enhanced features with/without automatic function	○
Castors	
Integral S single castors, (Ø 125 mm) – lying surface height = 31 – 86 cm	✓
Integral single castors, (Ø 150 mm) – lying surface height = 33.5 – 88.5 cm	○
Integral S single castors, (Ø 150 mm) – lying surface height = 33.5 – 88.5 cm	○
Linea twin castors, (Ø 125 mm) – lying surface height = 28.5 – 83.5 cm	○
Twin castors, (Ø 150 mm) – lying surface height = 34 – 90 cm	○
Linea twin castors, (Ø 150 mm) – lying surface height = 31 – 86 cm	○
Track caster	
at the foot of the bed	✓
at the head of the bed	○
5th castor, only in combination with twin castor (Ø 150 mm) lying surface height = 35 – 90 cm	○
Central brake	
Foot pedal at foot end, on both sides	✓
Foot pedal at foot and head end, on both sides	○
Additional configuration options	
Bed extension	○
Angle display – on both sides	○
Underbed light	○
Brake alarm	○
Attachment points on the lying surface for fixation belt holders	○

✓ Standard ○ optional